

# Determining Optimal Seeding Density for Human T Cell Expansion in G-Rex

David Hermanson<sup>2</sup>, Josh Ludwig<sup>1</sup>

1. ScaleReady, New Brighton, MN 55112  
2. Bio-Techne, Minneapolis, MN 55416

## Introduction


Efficiently expanding human T cells in a closed-system while maintaining a consistent cell phenotype is an active challenge for autologous and allogeneic T cell therapy manufacturers. Cell phenotype drift and variability in cell yield are common pain-points when transitioning from research protocols (i.e. flask or well plates) into larger scale, closed manufacturing systems (i.e. gas permeable bags, stirred tank bioreactors). This variability can be partially attributed to the introduction of motion and liquid media exchanges to accommodate cell expansion within larger volumes. G-Rex bioreactors have been shown (e.g. TIL, NK cells, T cells) to resolve these challenges by providing a static growth environment for cells, void of motion or media exchange, and promoting consistent cell expansion with reproducible cell phenotypes.

Optimizing cell expansion in G-Rex at small-scale is key to achieving efficient cell expansion. In this study we investigated the kinetics of human T cell expansion in small scale G-Rex6M bioreactors using ScaleReady reagents. Specifically, we determine the affect of seeding density on fold expansion and total cell yield. In addition, we look at the phenotypic composition of expanded T cells, including total CD4<sup>+</sup> and CD8<sup>+</sup> cells as well as memory CD4<sup>+</sup> T cells.

## Methods

Human CD3<sup>+</sup> T cells from 3 donors seeded at two densities (0.5x10<sup>6</sup>, 1.0x10<sup>6</sup>) into G-Rex<sup>®</sup>6M Well Plates. For each donor, separate G-Rex6M Well Plates were seeded for undisturbed cell expansion until harvest and cell analysis at Day 0, 3, 6, 8, and 10. Triplicates for each donor were run at each collection timepoint. All cells were cultured using ExCellerate<sup>™</sup> T Cell Expansion Media supplemented with GMP rhIL-7 and GMP rhIL-15. Cloudz<sup>™</sup> T Cell Activation Kit was used for cell activation. Cloudz CD3/CD28 micropsheres were added at Day 0 and removed at Day 2. After Cloudz removal, cell were returned to G-Rex 6M Well Plates for undisturbed rapid expansion using ExCellerate T Cell Expansion Media supplemented with IL-7 and IL-15. A detailed experimental procedure overview is depicted below.

### Experimental Overview



### Conditions

| Condition | Cell Number       | Cell Density                         | Cloudz Concentration          | Checkpoint / Harvest Day          | Glucose & Lactate Measurement              |
|-----------|-------------------|--------------------------------------|-------------------------------|-----------------------------------|--------------------------------------------|
| #1        | 5x10 <sup>6</sup> | 0.5x10 <sup>6</sup> /cm <sup>2</sup> | 25 µL/1x10 <sup>6</sup> cells | Day 3<br>Day 6<br>Day 8<br>Day 10 | Day 6<br>Day 7<br>Day 8<br>Day 9<br>Day 10 |
| #2        | 1x10 <sup>7</sup> | 1x10 <sup>6</sup> /cm <sup>2</sup>   | 25 µL/1x10 <sup>6</sup> cells |                                   |                                            |




\*Each condition was repeated with 3 donors

### Raw Materials


| Product Name                              | Concentration                 | Catalog Number | Supplier    |
|-------------------------------------------|-------------------------------|----------------|-------------|
| G-Rex <sup>®</sup> 6M Well Plate          | N/A                           | 80660M         | Wilson Wolf |
| GMP rhIL-7                                | 10 ng/mL                      | 207-GMP        | Bio-Techne  |
| GMP rhIL-15                               | 5 ng/mL                       | 247-GMP        | Bio-Techne  |
| Cloudz <sup>™</sup> T Cell Activation Kit | 25 µL/1x10 <sup>6</sup> cells | CLD001         | Bio-Techne  |
| ExCellerate <sup>™</sup> T Cell Media     | 100 mL/well                   | CCM030         | Bio-Techne  |

## Results

### Similar Cell Expansion Kinetics and Cell Yields Across Seeding Densities




### Higher Fold Expansion Using Lower Seeding Density




## Results

### High Fold Expansion of CD4<sup>+</sup> and CD8<sup>+</sup> Cells



### Rapid Expansion of Memory CD4<sup>+</sup> T Cell Population



## Conclusion

- Maximum T cell yield and expansion kinetics can be achieved using a lower starting cell density in the G-Rex. This improves the efficiency of T cell expansion by minimizing the number of donor cells needed to achieve maximum expansion within a G-Rex.
- G-Rex facilitates rapid expansion of T cells
- G-Rex T cell expansion is compatible with ScaleReady reagents:
  - Cloudz T Cell Activation Kit
  - GMP IL-7 & GMP IL-15
  - ExCellerate T Cell Expansion Media
- Small-scale preliminary data suggests:
  - Successful expansion of pre-selected CD4<sup>+</sup> and CD8<sup>+</sup> cells in a G-Rex6M Well Plate
  - Seeding densities of 0.5x10<sup>6</sup> cells/cm<sup>2</sup> and 1x10<sup>6</sup> cells/cm<sup>2</sup> produced 3.7x10<sup>8</sup> and 3.8x10<sup>8</sup> amount of total viable cells, respectively, with > 85% viability.
  - Desired yield, phenotype, and composition of final product may dictate harvest date.

Contact us: [info@scaleready.com](mailto:info@scaleready.com)

ScaleReady is a Joint Venture formed by Bio-Techne, Fresenius Kabi, and Wilson Wolf. Combining selected offerings from the three partners, the ScaleReady manufacturing platform combines tools and technologies for cell culture, cell activation and expansion, gene editing, and cell processing.

Trademarks and registered trademarks are the property of ScaleReady or parent companies unless otherwise specified.