

Generation of allogeneic CAR-NKT cells from hematopoietic stem and progenitor cells using a clinically guided culture method

Received: 23 August 2023

Accepted: 28 March 2024

Published online: 14 May 2024

Check for updates

Yan-Ruide Li^{1,18}, Yang Zhou^{1,18}, Jiaji Yu^{1,18}, Yu Jeong Kim¹, Miao Li¹, Derek Lee¹, Kuangyi Zhou¹, Yuning Chen¹, Yichen Zhu¹, Yu-Chen Wang², Zhe Li¹, Yanqi Yu¹, Zachary Spencer Dunn³, Wenbin Guo⁴, Xinjian Cen¹, Tiffany Husman¹, Aarushi Bajpai¹, Adam Kramer¹, Matthew Wilson¹, Ying Fang¹, Jie Huang¹, Shuo Li⁵, Yonggang Zhou⁵, Yuchong Zhang¹, Zoe Hahn¹, Enbo Zhu¹, Feiyang Ma⁶, Calvin Pan², Aldons J. Lusic^{1,2,7}, Jin J. Zhou⁸, Christopher S. Seet^{9,10,11,12}, Donald B. Kohn^{1,9,13}, Pin Wang³, Xianghong Jasmine Zhou⁵, Matteo Pellegrini^{4,14,15}, Benjamin R. Puliafito¹⁶, Sarah M. Larson^{10,17} & Lili Yang^{1,9,10,12}✉

Cancer immunotherapy with autologous chimeric antigen receptor (CAR) T cells faces challenges in manufacturing and patient selection that could be avoided by using ‘off-the-shelf’ products, such as allogeneic CAR natural killer T (AlloCAR-NKT) cells. Previously, we reported a system for differentiating human hematopoietic stem and progenitor cells into AlloCAR-NKT cells, but the use of three-dimensional culture and xenogeneic feeders precluded its clinical application. Here we describe a clinically guided method to differentiate and expand IL-15-enhanced AlloCAR-NKT cells with high yield and purity. We generated AlloCAR-NKT cells targeting seven cancers and, in a multiple myeloma model, demonstrated their antitumor efficacy, expansion and persistence. The cells also selectively depleted immunosuppressive cells in the tumor microenvironment and antagonized tumor immune evasion via triple targeting of CAR, TCR and NK receptors. They exhibited a stable hypoimmunogenic phenotype associated with epigenetic and signaling regulation and did not induce detectable graft versus host disease or cytokine release syndrome. These properties of AlloCAR-NKT cells support their potential for clinical translation.

CAR-T cell therapies¹ have received Food and Drug Administration approval for the treatment of B cell malignancies and multiple myeloma (MM)². However, autologous CAR-T cell products are associated with high costs, prolonged manufacturing time and limited accessibility to patients³. In particular, patients with progressive disease or who have received previous treatments may not have adequate or functional T cells for CAR-T cell production. To develop off-the-shelf cellular products, two approaches are being explored. One uses conventional $\alpha\beta$ T cells with abrogation of endogenous T cell receptor

(TCR) expression to minimize the risk of graft-versus-host disease (GvHD)⁴. The other uses cell types that inherently pose low GvHD risk, such as macrophages, NK cells and invariant natural killer T (iNKT or NKT) cells^{5–7}.

NKT cells, characterized by an invariant TCR α chain (TRAV10-TRAJ18 in human and TRAV11-TRAJ18 in mice) and coexpression of NK markers, are a unique type of $\alpha\beta$ T cell that is relatively rare in the human bloodstream, accounting for only 0.001–1% of circulating T cells^{8–11}. Owing to their TCR recognition of the non-polymorphic

A full list of affiliations appears at the end of the paper. ✉ e-mail: liliyang@ucla.edu

major histocompatibility complex-like molecule CD1d (ref. 12), these cells are not expected to induce GvHD¹³. NKT cells have demonstrated promising characteristics for the development of off-the-shelf cell therapy. They possess potent tumor-killing activity and can infiltrate tumors, and they bridge innate and adaptive immune responses¹⁴. In clinical studies, autologous CAR-engineered NKT (CAR-NKT) cells have shown efficacy in treating relapsed or resistant neuroblastoma without causing noticeable toxicity or cytokine release syndromes (CRSs)^{6,15}. However, the development of NKT cell therapies is limited by the scarcity of NKT cells in human blood and the difficulty of expanding them from peripheral blood mononuclear cells (PBMCs) while excluding conventional T cells, which could trigger GvHD¹³. Thus, identifying clinically adaptable methods to allogeneic NKT cells from sources other than PBMCs is of interest.

Hematopoietic stem and progenitor cells (HSPCs) have been genetically engineered and differentiated into various types of immune cell^{16,17}, including NKT cells^{18,19}. In a previous study, we differentiated engineered cord-blood HSPCs into NKT cells using an artificial thymic organoid culture method¹⁸. However, the use of three-dimensional (3D) culture and xenogeneic feeder cells precluded scale-up and clinical manufacturing¹⁸. To develop a clinically relevant system, here we describe a method for differentiating and expanding allogeneic CAR-NKT (AlloCAR-NKT) cells free of 3D culture and xenogeneic feeder cells. We generated AlloCAR-NKT cells targeting two blood cancer markers (BCMA and CD19) and three solid tumor markers (GD2, GPC3 and EGFRvIII). We also engineered AlloCAR-NKT cells with human IL-15, which is well recognized for its role in enhancing effector function, promoting memory formation and mitigating exhaustion in T, NK and NKT cells^{20–22}. Finally, we conducted comprehensive preclinical studies of AlloCAR-NKT cells to evaluate their manufacturing, pharmacology, efficacy, mechanism of action, pharmacokinetics and/or pharmacodynamics (PK–PD), safety and immunogenicity.

Results

Ex vivo generation of AlloCAR-NKT cells

Cryopreserved human cord-blood-derived CD34⁺ HSPCs were obtained from commercial vendors and used to generate NKT cells in a five-stage, 6-week process (Fig. 1a–c).

At stage 0, freeze–thawed CD34⁺ HSPCs were transduced with a lentivector codelivering a human iNKT TCR, a designated CAR, as well as additional genes, and cultured over 48 hours in a conventional X-VIVO 15-based HSPC medium (Fig. 1a,b,d and Supplementary Fig. 1a,b)^{18,23}. All lentivectors tested showed robust HSPC transduction (Supplementary Fig. 2a). The iNKT TCR genes were cloned from a healthy donor and validated for its effective recognition of CD1d and lipid antigens^{18,23}. CARs were studied that target various tumor antigens such as the B cell maturation antigen (BCMA), CD19, disialoganglioside GD2 (GD2), glypican-3 (GPC3) and epidermal growth factor receptor variant III (EGFRvIII), and contain different costimulatory domains including

CD28 and 4-1BB (Fig. 1d and Supplementary Fig. 1a)^{24–28}. An example of an additional gene (AG) is *IL-15*, which has been shown to augment the in vivo persistence and antitumor reactivity of NKT cells^{6,26,29,30}.

Gene-engineered HSPCs were then cultured over roughly 6 weeks to generate a designated AlloCAR-NKT cell product: stage 1 HSPC expansion (2 weeks), stage 2 NKT differentiation (1 week), stage 3 NKT deep differentiation (1 week) and stage 4 NKT expansion (2 weeks) (Fig. 1a,b,e and Supplementary Fig. 1c). We implemented the entire culture in a feeder-free and serum-free manner for compatibility with clinical translation³¹. We also demonstrated two alternative strategies that use human feeder cells, either α -galactosylceramide (α GC)-loaded healthy donor PBMCs or K562-based artificial antigen-presenting cells (aAPCs), during stage 4 NKT expansion, as these human feeder cells have been validated for clinical development (Fig. 1b and Supplementary Fig. 2b)^{15,31}.

The differentiation of AlloCAR-NKT cells followed a typical T cell lineage commitment and NKT cell developmental path, from CD4⁺ CD8[−] double-negative to CD4⁺ CD8⁺ double-positive and eventually to CD8⁺ single-positive or double-negative stage (Fig. 1e,g and Supplementary Fig. 1c)³². The final AlloCAR-NKT cells lacked a CD4⁺ single-positive population, which is typically present in endogenous human NKT cells (Fig. 1e and Supplementary Fig. 1c)³². Further examination of CD8⁺ AlloCAR-NKT cells showed a predominant expression of the CD8 α / α form (~80%) over the CD8 α / β form (~20%), resembling the CD8 α / α versus CD8 α / β expression pattern in the endogenous NKT cells (Supplementary Fig. 3)³². In general, CD8 α / α and double-negative human NKT cells are considered of similar functionality; these cells are pro-inflammatory and exhibit heightened cytotoxicity, making them desirable for cancer immunotherapy³². Incorporation of a CAR gene or CAR and IL-15 genes did not interfere with the differentiation of allogeneic NKT cells and yielded allogeneic NKT cells that universally coexpressed the CAR (and IL-15) alongside the iNKT TCR (Fig. 1f and Supplementary Fig. 1c).

The yield and purity of AlloCAR-NKT cells were high. A more than 10⁶-fold expansion was achieved from the input CD34⁺ HSPCs to the output mature AlloCAR-NKT cells (Fig. 1c). In our cultures, we routinely generated roughly 10¹⁰ mature AlloCAR-NKT cells from 10⁴ input HSPCs (Supplementary Fig. 2c–f). Flow cytometry and single-cell TCR sequencing analyses revealed uniform expression of the transgenic iNKT TCRs in AlloCAR-NKT cells (Fig. 1e,h and Supplementary Fig. 1c). Randomly recombined endogenous $\alpha\beta$ TCRs were undetected, likely attributed to allelic exclusion induced by the overexpression of the transgenic TCRs (Fig. 1e,h)³³.

The AlloCAR-NKT cell production process was robust, spanning 15 cord-blood donors and various cargo genes (Fig. 1c–e and Supplementary Fig. 1). These cargo genes included CARs with different designs, as well as additional genes such as *IL-15* (Fig. 1d and Supplementary Fig. 1). The variability in cord-blood donors and cargo genes had no discernible impact on the high yield and purity of the final product cells (Fig. 1c).


Fig. 1 | AlloNKT and AlloCAR-NKT cells are produced at high yield, purity and robustness. **a**, Schematics showing the generation of F^{Allo}(CAR)-NKT cells. AG, additional gene; CB, cord blood. **b**, Table showing the stages of F^{Allo}(CAR)-NKT cell culture. **c**, Yield of F^{Allo}(CAR)-NKT cells (AlloNKT and AlloBCAR-NKT, $n = 6$; Allo¹⁵BCAR-NKT and Allo¹⁵CAR19-NKT, $n = 9$; Allo¹⁵GD2.CAR-NKT, $n = 3$; Allo¹⁵GPC3.CAR-NKT and Allo¹⁵EGFRvIII.CAR-NKT, $n = 4$; n indicates different donors). **d**, Schematics showing the design of indicated lentivectors. Δ LTR, self-inactivating long-term repeats; MNDU3, internal promoter derived from the MND retroviral LTR U3 region; ϕ , packaging signal with the splicing donor and splicing acceptor sites; RRE, reverse-responsive element; cPPT, central polypurine tract; WPRE, woodchuck responsive element. **e**, FACS monitoring of the generation of F^{Allo}(CAR)-NKT cells. AlloBCAR-NKT, AlloNKT cells armed with BCAR; Allo¹⁵BCAR-NKT, AlloNKT cells armed with BCAR and IL-15. iNKT TCR was stained using a 6B11 monoclonal antibody. **f**, FACS detection of the iNKT TCR and BCAR on F^{Allo}(CAR)-NKT cells. BCAR was stained using an antimouse IgG F(ab')₂ antibody. **g**, UMAP visualization of

Allo¹⁵BCAR-NKT cell development. scRNA-seq data are presented. Each dot represents a single cell and is colored according to its cell cluster assignment. **h**, Single-cell TCR sequencing analyses of F^{Allo}(CAR)-NKT cells. Healthy donor PBMC-derived conventional $\alpha\beta$ T (PBMC-Tc), NKT (PBMC-NKT) and NK (PBMC-NK) cells were included as controls. **i**, Bulk RNA-seq analyses of F^{Allo}(CAR)-NKT cells. AlloNKT ($n = 3$), AlloBCAR-NKT ($n = 3$), Allo¹⁵BCAR-NKT ($n = 3$) and Allo¹⁵CAR19-NKT ($n = 2$) cells were analyzed. PBMC-Tc ($n = 8$), PBMC-NKT ($n = 3$) and PBMC-NK ($n = 2$) cells were included as controls. n indicates different donors. **j**, FACS detection of surface markers, intracellular cytokines and cytotoxic molecules in Allo¹⁵BCAR-NKT cells. PBMC-derived conventional T cells engineered with the same BCMA-targeting CAR (denoted as BCAR-T cells) were included as a control. Representative of one (**g–i**) and more than ten (**a–f** and **j**) experiments. Data are presented as the mean \pm s.e.m. and were analyzed by one-way ANOVA (**c**). NS, not significant.

In this study, we performed comprehensive preclinical characterization of $Allo^0$ CAR-NKT cells, with a primary focus on MM as the disease indication and its corresponding BCMA-targeting $Allo^0$ CAR-NKT cells, particularly those incorporating IL-15 enhancement ($Allo^{15}$ BCAR-NKT cells).

Phenotype and functionality of $Allo^0$ CAR-NKT cells

We analyzed the phenotype and functionality of $Allo^{15}$ BCAR-NKT cells in comparison with endogenous human PBMC-derived conventional $\alpha\beta$ T, NK and NKT cells (denoted PBMC-Tc, PBMC-NK and PBMC-NKT cells, respectively) (Supplementary Fig. 4a). PBMC-Tc cells engineered with



the same BCMA-targeting CAR (BCAR) (conventional BCAR-T cells) were also included as a benchmark control (Supplementary Fig. 4b,c)²⁴.

Gene profiling using bulk RNA sequencing (RNA-seq) revealed that $Allo^{15}$ BCAR-NKT cells closely resembled PBMC-NKT cells, ranked second to PBMC-Tc cells and were largely different from PBMC-NK cells (Fig. 1i). To validate iNKT TCR function, we stimulated $Allo^{15}$ BCAR-NKT cells with an agonist glycolipid antigen α GC (Supplementary Fig. 5). $Allo^{15}$ BCAR-NKT cells proliferated robustly (Supplementary Fig. 5a–c) and secreted high levels of Th1 cytokines (that is, IFN γ and TNF α) while low levels of Th2/Th17 cytokines (that is, IL-4 and IL-17A) (Supplementary Fig. 5d), indicating a Th1-prone functionality of $Allo^{15}$ BCAR-NKT cells, which agrees with their CD8 single-positive and/or double-negative phenotype regardless of IL-15 engineering (Fig. 1e and Supplementary Fig. 5e)³².

Flow cytometry analysis of $Allo^{15}$ BCAR-NKT cells revealed a typical NKT cell phenotype (Fig. 1j). Compared with conventional BCAR-T cells, $Allo^{15}$ BCAR-NKT cells expressed high levels of NK receptors (NKR) (for example, CD161) and central memory markers (for example, CD62L, CD134 and LEF1), and produced exceedingly high levels of effector cytokines (for example, IFN γ , IL-2 and TNF α) and cytotoxic molecules (for example, Granzyme B and Perforin), agreeing with their CD8 single-positive and/or double-negative phenotype (Fig. 1j and Supplementary Fig. 5f)^{34,35}.

In vitro antitumor efficacy and mechanism of action of $Allo^{15}$ BCAR-NKT cells

$Allo^{15}$ BCAR-NKT cells are expected to target tumor cells through multiple surface receptors, involving CAR recognition of CAR antigen, TCR recognition of CD1d and NKR recognition of NK ligands expressed on most tumor cells (Fig. 2a)¹⁸. $Allo^{15}$ BCAR-NKT cells could deploy CAR/TCR/NKR triple-targeting mechanisms, suitable for CD1d⁺ tumors, such as certain types of blood cancer (for example, MM, acute myeloid leukemia and myelomonocytic leukemia) and some solid tumors (for example, medulloblastoma and glioblastoma) (Fig. 2a)¹⁸. For CD1d⁻ tumors, including other blood cancers and most solid tumors, $Allo^{15}$ BCAR-NKT cells could still deploy CAR/NKR double-targeting mechanisms (Fig. 2a)¹⁸.

A cohort of primary bone marrow (BM) samples from patients with MM was collected and used to evaluate the tumor cell killing efficacy of $Allo^{15}$ BCAR-NKT cells (Fig. 2b). Coexpression of BCMA, CD1d and NK ligands was detected on MM cells from all samples (Fig. 2c and Supplementary Fig. 6). In an in vitro tumor cell killing assay, $Allo^{15}$ BCAR-NKT cells demonstrated superior effectiveness in eliminating primary MM cells compared with conventional BCAR-T cells (Fig. 2d,e). Even without CAR engineering, $Allo^{15}$ NKT cells displayed considerable tumor-killing capabilities, albeit slightly lower than that of $Allo^{15}$ BCAR-NKT cells (Fig. 2d,e), indicating the CAR-independent tumor-killing potential of $Allo^{15}$ BCAR-NKT cells.

The CAR/TCR/NKR triple-targeting mechanisms of $Allo^{15}$ BCAR-NKT cells were further validated using a series of in vitro tumor cell killing assays (Fig. 2f–n), involving three human

MM.1S-derived cell lines as targets (MM-FG-CD1d that was BCMA⁺CD1d⁺, MM-FG that was BCMA⁺CD1d⁻ and KO MM-FG that was BCMA⁻CD1d⁻) and four therapeutic cells ($Allo^{15}$ BCAR-NKT cells, and three controls cells: PBMC-T, BCAR-T and $Allo^{15}$ NKT) (Fig. 2f–h and Supplementary Fig. 7a–c). BCAR-T cells solely relied on BCAR-BCMA recognition for tumor cell cytotoxicity, as evidenced by selective killing of BCMA⁺ tumor cells (Fig. 2i). In contrast, $Allo^{15}$ BCAR-NKT cells exhibited cytotoxicity against both BCMA⁺ and BCMA⁻ tumor cells, indicating a less pronounced dependence on BCAR-BCMA recognition, albeit with a comparatively reduced killing capacity toward BCMA⁻ tumor cells (Fig. 2i). Meanwhile, $Allo^{15}$ BCAR-NKT cells killed tumor cells more effectively in the presence of CD1d/ α GC, indicating a TCR-directed targeting mechanism (Fig. 2i). Moreover, $Allo^{15}$ BCAR-NKT cells could kill BCMA⁻CD1d⁻ tumor cells through NKR (that is, NKG2D and DNAM-1) recognition, confirming an NKR-mediated targeting mechanism (Fig. 2k,l and Supplementary Fig. 7d–f).

We visualized the direct recognition and attack of MM.1S tumor cells by $Allo^{15}$ BCAR-NKT cells by scanning electron microscopy (SEM) (Fig. 2j). Corresponding to their multiple targeting mechanisms (Fig. 2f–l) and strong cytotoxic functions (Fig. 1j), $Allo^{15}$ BCAR-NKT cells displayed superior tumor cell killing efficacy compared with BCAR-T cells against both primary patient samples (Fig. 2e) and tumor cell lines (Fig. 2i), which was associated with upregulation of activation markers (that is, CD69) and production of cytotoxic molecules (that is, Perforin and Granzyme B) and effector cytokines (that is, IFN γ) (Fig. 2m,n and Supplementary Fig. 7g,h).

In vivo antitumor efficacy and PK–PD of $Allo^{15}$ BCAR-NKT cells

To investigate the in vivo antitumor efficacy and PK–PD of $Allo^{15}$ BCAR-NKT cells, we conducted a pair of mirror-imaging experiments using a human MM xenograft NSG (NOD scid gamma) mouse model (Fig. 3)¹⁸. Conventional BCAR-T cells were included as a benchmark control. In the efficacy experiment, MM tumor cells were labeled with firefly luciferase and enhanced green fluorescence protein (FG) dual reporters (Fig. 3a–e), while in the PK–PD experiment, therapeutic cells were labeled with the FG dual reporters (Fig. 3f–k and Supplementary Fig. 8a,b). This experimental setup enabled simultaneous visualization of the in vivo dynamics and biodistribution of both tumor cells and therapeutic cells during the same antitumor response.

Under a heavy tumor load condition¹⁸, a single administration of $Allo^{15}$ BCAR-NKT cells achieved tumor elimination in most (nine out of ten) experimental mice and long-term (>80 days) survival of all experimental mice (Fig. 3c–e). In contrast, both BCAR-T and $Allo^{15}$ BCAR-NKT cells only achieved partial suppression of tumor growth and limited improvement of survival (Fig. 3c–e). The experimental mice treated with BCAR-T cells died of a combination of tumor growth and GvHD (Fig. 3e)¹⁸.


In the PK–PD study under the same heavy tumor load condition, $Allo^{15}$ BCAR-NKT/FG cells expanded more than 200-fold, peaking around week 2, followed by a gradual decline and a long-term persistence

Fig. 2 $Allo^{15}$ BCAR-NKT cells directly kill tumor cells and use multiple targeting mechanisms. **a**, Diagram showing the mechanisms used by $Allo^{15}$ BCAR-NKT cells to kill CD1d⁺ or CD1d⁻ tumor cells. AML, acute myeloid leukemia; JMML, juvenile myelomonocytic leukemia; ALL, acute lymphoblastic leukemia; HL, Hodgkin lymphoma; NHL, non-Hodgkin lymphoma; MB, medulloblastoma; GBM, glioblastoma. **b–e**, Studying the antitumor efficacy of $Allo^{15}$ BCAR-NKT cells against primary samples from patients with MM. **b**, Diagram showing the collection of BM samples from patients with MM. **c**, FACS detection of CAR target (BCMA), iNKT TCR target (CD1d) and NKR target (ULBP-1) on primary patient-derived tumor cells. **d**, Experimental design to study the primary MM tumor cell killing by therapeutic cells. **e**, Tumor cell killing data at 24 h ($n = 3$). **f–j**, Studying the antitumor efficacy of $Allo^{15}$ BCAR-NKT cells against human MM.1S cell lines. **f**, Experimental design. **g**, Schematics showing the indicated human MM.1S cell lines. MM-FG, MM.1S cell line engineered to express the firefly luciferase and

enhanced green fluorescence protein (FG) dual reporters; MM-FG-CD1d, MM-FG cell line further engineered to overexpress human CD1d; KO MM-FG, MM-FG cell line further engineered to knockout the BCMA encoding gene. **h**, Schematics showing the indicated therapeutic cells. **i**, Tumor cell killing data at 24 h ($n = 4$). **j**, SEM image showing that an $Allo^{15}$ BCAR-NKT cell (blue) is attacking a MM.1S-FG tumor cell (purple). Scale bar, 5 μ m. **k–l**, Studying the tumor-killing mechanism of $Allo^{15}$ BCAR-NKT cells mediated by NK activating receptors (that is, NKG2D and DNAM-1). **k**, Experimental design. **l**, Tumor cell killing data at 24 h (E:T ratio 10:1; $n = 4$ from four different cell product donors). **m,n**, FACS characterization of $Allo^{15}$ BCAR-NKT cells 24 h after coculturing with MM-FG. **m**, FACS detection of surface CD69 as well as intracellular Perforin and Granzyme B in $Allo^{15}$ BCAR-NKT cells. **n**, Quantification of **m** ($n = 3$ from three different cell product donors). Representative of three experiments. Data are presented as the mean \pm s.e.m. and were analyzed by one-way ANOVA (**e**, **l** and **n**).

greater than 80 days (Fig. 3h,i). Additionally, these cells also exhibited strong tumor homing capacity, effectively and preferentially homing to the major MM tumor site, the BM, during both their peak response and extended persistence (Fig. 3j,k). This homing capacity likely contributes to their superior response rate (Fig. 3c–e). In the absence of IL-15 enhancement, the ^{Allo}BCAR-NKT/FG cells showed a moderate fivefold expansion, peaking around week 1, followed by a rapid decline and a greatly shortened persistence of around 10 days (Fig. 3h,i,k).

Conventional BCAR-T/FG cells displayed a moderate expansion and a rapid decline similar to those of ^{Allo}BCAR-NKT/FG cells in the first 40 days (Fig. 3h,i), which is likely attributed to the stimulation driven by tumor antigens². However, after 40 days (about 1 month after BCAR-T/FG adoptive transfer), BCAR-T/FG cells underwent a swift expansion and infiltrated into various major organs, leading to serious GvHD and mouse demise (Fig. 3h–k), consistent with observations in many CAR-T cell xenograft models^{18,36,37}.






Fig. 3 | **AlloCAR-NKT cells exhibit potent antitumor efficacy in vivo, associated with effective tumor homing, expansion and long-term persistence.**
a–e, Studying the in vivo antitumor efficacy of $Allo^{15}$ BCAR-NKT cells in a human MM xenograft NSG mouse model. MM-FG indicates MM.1S tumor cell line labeled with the FG dual reporters. **a**, Experimental design. **b**, Schematics showing the tumor cell line and the therapeutic cells used in the study. **c**, BLI images showing the presence of tumor cells in experimental mice over time. **d**, Quantification of **c** (Vehicle, $n = 9$; BCAR-T and $Allo^{15}$ BCAR-NKT, $n = 10$). TBL, total body luminescence. p/s, photons per second. **e**, Kaplan–Meier survival curves of experimental mice over time (Vehicle, $n = 9$; BCAR-T and $Allo^{15}$ BCAR-NKT, $n = 10$). **f–k**, Studying the in vivo PK–PD of $Allo^{15}$ BCAR-NKT cells in a human MM xenograft NSG mouse model.

MM indicates MM.1S tumor cell line. Note the therapeutic cells, but not the MM.1S tumor cell line, were labeled with the FG dual reporters. **f**, Experimental design. **g**, Schematics showing the tumor cell lines and the therapeutic cells used in the PK–PD study. **h**, BLI images showing the presence of therapeutic cells in experimental mice over time. **i**, Quantification of **h** (BCAR-T/FG, $n = 3$; $Allo^{15}$ BCAR-NKT/FG, $n = 4$). **j**, BLI images showing the biodistribution of $Allo^{15}$ BCAR-NKT cells in a representative experimental mouse. Ventral, left and dorsal views are shown. **k**, Quantification of therapeutic cell tissue distribution (BCAR-T/FG, $n = 3$; $Allo^{15}$ BCAR-NKT/FG, $n = 4$). Representative of three experiments. Data are presented as the mean \pm s.e.m. and were analyzed by one-way ANOVA (**k**). min., minimum; max., maximum; ND, non-detectable; FD, found dead.

In vivo gene profile of $Allo^{CAR}$ -NKT cells

To study the genomic and molecular regulation of the in vivo antitumor performance of $Allo^{15}$ BCAR-NKT cells, we used a human MM xenograft NSG mouse model in combination with single-cell RNA sequencing (scRNA-seq) analyses (Fig. 4a). Conventional BCAR-T cells were included as a benchmark control. Three types of therapeutic cell ($Allo^{15}$ BCAR-NKT, $Allo^{15}$ BCAR-NKT and BCAR-T cells) were collected at three time points (days 0, 28 and 35) followed by scRNA-seq analysis (Fig. 4a). Day 0 samples were collected before infusion, representing the starting status of the therapeutic cells. Day 28 samples were collected from the BM tumor sites of the experimental mice, representing the in vivo tumor-responding status of the therapeutic cells. Day 35 samples were generated by stimulating day 28 samples in vitro with BCMA-expressing aAPCs for 1 week, representing the reactivation potential of the therapeutic cells posttumor exposure (Fig. 4a).

Uniform manifold approximation and projection (UMAP) analysis of combined samples revealed the formation of five major cell clusters (Fig. 4b). The signature gene profiling and gene set enrichment analysis identified cluster 1 cells as proliferating cells, cluster 2 cells as effector and/or memory cells, cluster 3 cells as early exhausted cells, cluster 4 cells as terminally exhausted cells and cluster 5 cells as resting cells (Fig. 4b and Supplementary Fig. 9a,b)^{38–45}.

All three types of therapeutic cell displayed large gene profile changes when comparing their day 0 preinfusion samples with their day 28 in vivo samples, suggesting an in vivo genomic and/or molecular signature of a cell product that cannot be fully captured by analyzing its preinfusion product (Fig. 4c,d), consistent with multiple preclinical and clinical studies^{46,47}. Day 35 restimulated $Allo^{15}$ BCAR-NKT cells exhibited a tumor cell killing capacity similar to their day 0 preinfusion counterparts, in contrast to the reduced tumor cell killing capability observed in day 35 restimulated BCAR-T cells (Fig. 2i left panel and Supplementary Fig. 10a,b). This disparity aligns with a better-preserved effector and/or memory cell population (cluster 2) observed in restimulated $Allo^{15}$ BCAR-NKT cells, compared with that observed in the restimulated BCAR-T cells (Fig. 4c,d). These results indicate that allogeneic CAR-NKT cells, despite their prolonged in vivo exposure to tumors, retain substantial potential for reactivation that may result in sustained therapeutic benefits⁴⁸.

We conducted detailed analyses of the day 28 in vivo samples. Our PK–PD study indicated that day 28 falls in the declining phase of $Allo^{15}$ BCAR-NKT cell expansion and in the initiation phase of the second BCAR-T cell expansion, likely due to GvHD (Fig. 3h,i). Accordingly, BCAR-T cells contained a notable population of proliferating cells (cluster 1), which was less prevalent in $Allo^{15}$ BCAR-NKT cells (Fig. 4c,d). When we focused on analyzing the non-proliferating cell populations (clusters 2, 3 and 4), compared with BCAR-T and $Allo^{15}$ BCAR-NKT cells, $Allo^{15}$ BCAR-NKT cells contained more effector and memory cells (cluster 2), fewer early exhausted cells (cluster 3) and a similar fraction of terminally exhausted cells (cluster 4) (Fig. 4e). In support of this notion, RNA-velocity analysis revealed the accumulation of $Allo^{15}$ BCAR-NKT cells

at the effector and/or memory stage, while BCAR-T and particularly $Allo^{15}$ BCAR-NKT cells showed retention at the early exhausted stage (Fig. 4f).

Pathway and gene signature analyses showed that compared with BCAR-T cells, $Allo^{15}$ BCAR-NKT cells displayed enhanced effector and memory, and attenuated exhaustion gene profiles (Fig. 4g,i), aligning with their superior in vivo antitumor efficacy (Fig. 3c–e). These features are likely attributed to IL-15 expression, as they were absent from non-IL-15-engineered $Allo^{15}$ BCAR-NKT cells (Figs. 3c–e and 4g,i)^{15,22,30,49}. Consequently, an augmented IL-15 signaling pathway was observed in $Allo^{15}$ BCAR-NKT cells but not in BCAR-T or $Allo^{15}$ BCAR-NKT cells (Fig. 4g).

Nevertheless, there are fundamental differences between the allogeneic CAR-NKT and conventional CAR-T cell platforms. In comparison to BCAR-T cells, both $Allo^{15}$ BCAR-NKT and $Allo^{15}$ BCAR-NKT cells exhibited fewer T cell features and markedly elevated NK cell features, particularly in the in vivo day 28 samples (Fig. 4h,i). The NK cell features of $Allo^{15}$ BCAR-NKT cells were reflected in their heightened expression of NKR and robust NKR-mediated tumor cell killing (Fig. 2k,l and Supplementary Fig. 7d,f). BCAR-T cells acquired certain NK cell feature in vivo (Fig. 4h), consistent with a previous study reporting an NK-like transition and dysfunction in CAR-T cells subjected to repetitive tumor stimulation⁴⁶. Conversely, in vivo augmentation of NK features in $Allo^{15}$ BCAR-NKT cells did not correlate with signs of exhaustion (Fig. 4e,g,i).

Antagonize tumor immune evasion by $Allo^{CAR}$ -NKT cells

In addition to analyzing the in vivo gene profiles of therapeutic cells, we conducted scRNA-seq analyses of the MM tumor cells collected from the same in vivo experiment (Figs. 4a and 5a). Our primary focus was to investigate tumor immune evasion, in particular the loss of CAR antigen, which is considered an obstacle for CAR-directed cell therapy⁵⁰.

We observed a reduction of BCMA antigen expression in tumors treated with BCAR-T cells (Fig. 5b,c and Supplementary Fig. 11), a phenomenon that has been reported in patients with MM under BCAR-T cell therapy⁵¹. In contrast, no BCMA antigen loss was observed in tumors treated with $Allo^{15}$ BCAR-NKT cells (Fig. 5b,c and Supplementary Fig. 11), which could be attributed to the multiple tumor-targeting mechanisms employed by $Allo^{15}$ BCAR-NKT cells (Fig. 2).

$Allo^{CAR}$ -NKT cells target the TME


The immunosuppressive tumor microenvironment (TME) is considered a notable barrier hindering cancer immunotherapy, particularly in the context of solid tumors, but also for certain blood cancers such as MM, especially when the disease involves BM sites^{52–54}. Tumor-associated macrophages and myeloid-derived suppressor cells (MDSCs), two major components of the immunosuppressive TME⁵⁴, both express high levels of CD1d, making them direct targets of NKT cells⁵⁴. We studied the interplay between $Allo^{15}$ BCAR-NKT cells and the MM TME using primary samples from patients with MM (Fig. 5d) and an in vivo MM NSG xenograft model (Fig. 5h).

Fig. 4 | $Allo^{CAR}$ -NKT cells display an in vivo gene profile associated with mixed T/NK cell features, strong effector and/or memory function and attenuated exhaustion property. **a**, Schematic showing the experimental design to study the gene profiling of $Allo^{15}$ BCAR-NKT cells using scRNA-seq. Conventional BCAR-T cells and non-IL-15 engineered $Allo^{15}$ BCAR-NKT were included as two controls. Therapeutic cell samples were collected at three time points: day 0 (preinfusion samples), day 28 (in vivo samples) and day 35 (in vitro restimulated samples), followed by FACS and scRNA-seq analyses. **b**, Combined UMAP plot showing the formation of five major cell clusters. Total cells combined from all samples are included. Each dot represents a single cell and is colored according to its cell cluster assignment. **c**, Individual UMAP plots showing cell cluster composition of the indicated samples. Each dot represents a single cell and is colored according to its cell cluster assignment. **d**, Bar graphs showing the cell cluster proportions of the indicated samples. **e**, Pie charts showing the cell cluster proportions of the

day 28 samples, focusing on cell clusters 2, 3 and 4. **f**, RNA-velocity projected on UMAP plots. Arrows represent the estimates of local average velocity, showing the path (indicated by arrow orientation) and pace (indicated by arrow length) of cell transition. Diagrams summarizing the UMAP cell clustering and RNA-velocity analysis are also presented. The cell cluster proportions are represented by the sizes of circles, while the transitions between cell clusters are symbolized using connecting arrows. **g**, Pathway analyses of differentially expressed genes comparing $Allo^{15}$ BCAR-NKT cells with BCAR-T cells (top) or $Allo^{15}$ BCAR-NKT cells (bottom). **h**, Violin plots showing the expression distribution of T and NK cell signature genes in the indicated therapeutic cell samples. **i**, Heatmap showing the expression of representative signature genes. Each column indicates an individual cell. Each row indicates an individual gene. The experiment was performed once; cells isolated from ten mice of each experimental group were combined for analyses.

In the first study, we cocultured $Allo^{15}$ BCAR-NKT cells with primary BM samples from patients with MM (Fig. 5d,e). $Allo^{15}$ BCAR-NKT cells effectively and selectively depleted tumor-associated macrophages and MDSCs, corresponding to the high CD1d expression on these

cells (Fig. 5f,g). $Allo^{15}$ BCAR-NKT cells spared other immune cell populations expressing no or low levels of CD1d, including granulocytes, T cells, B cells, NK cells and HSPCs (Fig. 5f,g). Further studies focused on the subpopulations of HSPCs isolated from either BM samples



from patients with MM or G-CSF-mobilized healthy donor leukopaks, including long-term HSPCs, short-term HSPCs and multi-potent progenitor cells, showing that all three HSPC subpopulations expressed undetectable CD1d and were not subjected to $Allo^{15}$ BCAR-NKT cell killing (Supplementary Fig. 12a–j).

In the second study, we administered $Allo^{15}$ BCAR-NKT cells into NSG mice bearing human MM xenografts, and collected BM cells and other tissues after 28 days for analysis (Fig. 5h). BM cells collected from tumor-free NSG mice (blank NSG), non-treated tumor-bearing NSG mice (NT) and BCAR-T cell-treated tumor-bearing NSG mice (BCAR-T) were included as controls (Fig. 5h). In this xenograft mouse model, the presence of cross-species reactivity between the human iNKT TCR and mouse CD1d allowed us to study the $Allo^{15}$ BCAR-NKT cell in vivo targeting of mouse CD1d-positive cells within the TME as well as across other tissues⁵⁵. scRNA-seq analysis of BM samples revealed the formation of seven major mouse cell clusters, encompassing HSPCs, dendritic cells (DCs), monocytes, M2-like macrophages, CRS-related macrophages, MDSCs and granulocytes (Fig. 5i and Supplementary Fig. 13). Engraftment of human MM cells led to an increase in the proportion of TME signature cells, such as the M2-like macrophages, CRS-related macrophages and MDSCs (Fig. 5j,k)^{54,56–58}. $Allo^{15}$ BCAR-NKT cell treatment altered the TME cell composition, selectively depleting TME signature cells while sparing other cells (that is, HSPCs, DCs, monocytes and granulocytes) (Fig. 5j–n). All three mouse HSPC subpopulations were preserved, including long- and short-term HSPCs and multi-potent progenitor cells (Supplementary Fig. 12k–m). In contrast, the BCAR-T cell treatment exacerbated a more hostile TME, reflected by a further enrichment of TME signature cells (Fig. 5j–n).

The selective $Allo^{15}$ BCAR-NKT cell targeting of BM TME signature cells (M2-like macrophages, CRS-related macrophages and MDSCs) was correlated with their high expression of CD1d, as evidenced by flow cytometry analyses (Fig. 5o). $Allo^{15}$ BCAR-NKT cells did not deplete residential monocytes and/or macrophages and granulocytes from other tissues (that is, peripheral blood, liver and peritoneal fluid), including the CD1d^{hi} tissue-resident macrophages found in the peritoneal fluid (Fig. 5o). This feature may arise from the preferential localization of $Allo^{15}$ BCAR-NKT cells to tumor sites (Fig. 3h,j,k), along with the absence of pronounced CD1d expression on quiescent monocytes, macrophages and/or granulocytes⁵².

Safety of $Allo^{15}$ BCAR-NKT cells

CRS-related macrophages have been reported to exacerbate CRS^{59,60}. Our scRNA-seq study identified a CRS-related macrophage population (cluster 5) in the BM of MM-engrafted mice (Fig. 5i). Treatment with $Allo^{15}$ BCAR-NKT cells, but not BCAR-T cells, effectively depleted this CRS-related macrophage population in the BM (Fig. 5j,k,n) and reduced CRS-related biomarker (that is, IL-6 and SAA-3) measurements in the serum (Fig. 5p,q)^{59,60}.

To evaluate GvHD risk, we used a traditional in vitro mixed lymphocyte reaction (MLR) assay (Fig. 6a) as well as an in vivo MM xenograft model (Fig. 3a). In the in vitro MLR assay, $Allo^{15}$ BCAR-NKT cells were stimulated with irradiated PBMCs from 15 mismatched healthy donors (Fig. 6a,b). Minimal IFN γ production was detected, whereas BCAR-T cells vigorously produced IFN γ (Fig. 6b).

In the in vivo MM xenograft model, BCAR-T cells elicited severe GvHD, leading to the death of experimental animals associated with the pronounced immune cell infiltration of vital organs such as the heart, liver and lung (Figs. 3e and 6c,d). In contrast, $Allo^{15}$ BCAR-NKT cell treatment achieved GvHD-free, long-term survival of experimental animals, with no signs of immune cell infiltration in vital organs (Figs. 3e and 6c,d).

Hypoimmunogenicity of $Allo^{15}$ BCAR-NKT cells


To evaluate allorejection of $Allo^{15}$ BCAR-NKT cells, we conducted two in vitro MLR assays (Fig. 6c,g). The first MLR assay was designed to study T cell-mediated allorejection, wherein irradiated $Allo^{15}$ BCAR-NKT cells were cocultured with PBMCs from 15 mismatched healthy donors, followed by measuring IFN γ production (Fig. 6e). Compared with BCAR-T cells, $Allo^{15}$ BCAR-NKT cells induced lower production of IFN γ (Fig. 6f). The second MLR assay was designed to study NK cell-mediated allorejection, wherein $Allo^{15}$ BCAR-NKT cells were cocultured with PBMC-NK cells from ten mismatched healthy donors, followed by quantifying viable $Allo^{15}$ BCAR-NKT cells (Fig. 6g). Compared with BCAR-T cells, $Allo^{15}$ BCAR-NKT cells showed greatly improved survival (Fig. 6h). Notably, while HSPC-engineered NKT cells were similar to endogenous PBMC-NKT cells in terms of displaying low GvHD risk, their distinctive characteristic of resistance to allorejection is unique, because the endogenous PBMC-NKT cells induced high production of IFN γ by donor-mismatched PBMCs in the first MLR assay, and were killed by donor-mismatched PBMC-NK cells in the second MLR assay (Supplementary Fig. 14).

To investigate the mechanisms responsible for the unique allorejection resistance of $Allo^{15}$ BCAR-NKT cells, we started by examining cell surface molecules associated with allorejection¹⁸. Compared with BCAR-T cells, $Allo^{15}$ BCAR-NKT cells expressed lower levels of surface HLA-I, HLA-II and NK ligands (that is, ULBP-1) in both in vitro cell cultures and in vivo antitumor studies (Fig. 6i–k). RNA-seq analyses revealed a sustained low expression of gene sets encoding HLA-I/II molecules and NK ligands (Fig. 6l,m), indicating an intrinsic and stable hypoimmunogenic phenotype of $Allo^{15}$ BCAR-NKT cells that is regulated at the transcriptional level. Next, we studied epigenetic regulation and signaling networks that govern transcription of the pertinent genes.

Chromatin remodeling processes regulate transcription of HLA and NK ligand genes^{61,62}. Genome-wide methylation sequencing (methyl-seq) analysis of $Allo^{15}$ BCAR-NKT cells revealed hypermethylation in the promoter regions of many HLA and NK ligand genes (Fig. 6n),

Fig. 5 | $Allo^{15}$ BCAR-NKT cells block tumor antigen escape, alter TME and exhibit low CRS characteristics. **a–c**, Studying the in vivo gene profile of MM tumor cells using scRNA-seq. **a**, Experimental design. **b**, UMAP plots showing the expression distribution of *TNFRSF17* gene encoding BCMA in the indicated samples. NT, no therapeutic cell treatment. **c**, Violin plot showing *TNFRSF17* gene expression. **d–g**, Studying the TME targeting using primary samples from patients with MM. **d**, Schematics showing the collection of primary samples from patients with MM. **e**, FACS analysis of immune cell composition. Gran, granulocyte; Mono, monocytes; Lympho, lymphocytes. **f,g**, In vitro culturing of samples from patients with MM with or without $Allo^{15}$ BCAR-NKT cells for 24 h. NT, no $Allo^{15}$ BCAR-NKT cell coculture. **f**, Killing data of the indicated TME cell component ($n = 3$). **g**, FACS measurement of surface CD1d expression in the indicated TME cell component. **h–o**, Studying the TME targeting using an in vivo human MM xenograft mouse model. **h**, Experimental design. **i**, Combined UMAP plot showing the formation of seven major cell clusters. DC, dendritic cell. **j**, Individual UMAP plots showing cell clusters of the indicated BM samples.

k, Bar graphs showing the cell cluster proportions of the indicated samples. **l–n**, Expression (expr.) of *Tgfb1* (**l**), M2-like M ϕ gene signature (**m**) and CRS-related M ϕ gene signature (**n**) in the indicated samples. **o**, FACS analyses of CD1d expression on mouse granulocytes and monocytes/macrophages. TRM, tissue-resident macrophage. **p,q**, Studying CRS response using an in vivo human MM xenograft mouse model. **p**, Experimental design. **q**, ELISA analyses of mouse IL-6 and SAA-3 in mouse serum ($n = 4$). SAA-3, Serum amyloid A-3. Representative of one (**a–c** and **h–n**) and three (**d–g** and **o–q**) experiments. For scRNA-seq, cells isolated from ten mice were combined for analyses. Data are presented as the mean \pm s.e.m. and were analyzed by two-tailed Student's *t*-test (**f**) or one-way ANOVA (**q**). *P* values of violin plots were determined by two-tailed Wilcoxon rank sum test (**l–n**) or Kruskal–Wallis test for the overall comparison and Dunn's test for post hoc pairwise comparisons between groups (**c**). In the violin plots, box and whisker plots exhibit the minimum, lower quartile, median, upper quartile and maximum expression levels of each type of cell.



suggesting that epigenetic regulation may contribute to the hypoinmunogenic phenotype of $Allo^{15}$ BCAR-NKT cells.

We examined the IFN γ -JAK-STAT1 signaling pathway, one of the best-characterized pathways regulating HLA and NK ligand gene expression^{63,64}. Compared with PBMC-T and BCAR-T cells, $Allo^{15}$ BCAR-NKT cells exhibited a desensitized IFN γ -JAK-STAT1 signaling, characterized by low STAT1 phosphorylation at both basal level and after IFN γ stimulation (Fig. 6o). Correspondingly, transcription factors encoded by IFN γ primary responsive genes (for example, *IRF1*) and secondary responsive genes (for example, *CIITA* and *SP-1*) were also maintained at low levels (Fig. 6o). IRF1 regulates *HLA-I* gene expression⁶⁵, CIITA regulates *HLA-II* expression⁶⁶ and SP-1 regulates NK ligands (for example, ULBP)⁶⁷. Overexpression of a constitutively active STAT1 molecule⁶⁸ partially reversed the hypoinmunogenic phenotype of $Allo^{15}$ BCAR-NKT cells, evidenced by their upregulation of surface HLA-I/II and NK ligands (Fig. 6p).

Collectively, these findings reveal a hypoinmunogenic phenotype that is regulated by epigenetic and signaling network mechanisms (Fig. 6q). This hypoinmunogenic phenotype is intrinsic to our allogeneic HSPC-engineered NKT cells, and is stable throughout inflammation (for example, IFN γ stimulation) and in vivo antitumor responses.

Comparison of HSPC- and PBMC-derived CAR-NKT cells

We performed side-by-side comparison of $Allo^{15}$ BCAR-NKT cells and PBMC-derived IL-15-enhanced BCAR-NKT ($PBMC^{15}$ CAR-NKT) cells (Figs. 7a–j). Flow cytometry analysis revealed substantial similarities between the two types of BCAR-NKT cell, including their capacity to produce high levels of effector cytokines (that is, IFN γ , TNF α and IL-2) and cytotoxic molecules (that is, Granzyme B and Perforin) (Fig. 7e). However, compared with $PBMC^{15}$ CAR-NKT cells, $Allo^{15}$ BCAR-NKT cells expressed higher levels of NKR (for example, DNAM-1, NKp44 and NKp30), indicating their enhanced NK feature (Fig. 7d). The antitumor efficacy of both BCAR-NKT cells was assessed using a series of in vitro tumor cell killing assays, as well as a human MM xenograft NSG mouse model (Fig. 7f,h). Both BCAR-NKT cells demonstrated similar antitumor efficacy in killing BCMA⁺ MM tumor cells, both in vitro and in vivo (Fig. 7g–k). Compared with $PBMC^{15}$ CAR-NKT cells, $Allo^{15}$ BCAR-NKT cells exhibited enhanced killing of BCMA⁺ tumor cells, particularly those sensitive to NKR-mediated cytotoxicity (for example, human melanoma cell line A375 and leukemia line K562), aligning with their heightened NKR expressions (Fig. 7d,g).

We integrated the scRNA-seq datasets reported from a recent clinical trial of autologous PBMC-derived GD2.CAR-engineered IL-15-enhanced NKT cells¹⁵ with those generated from our preclinical animal study evaluating $Allo^{15}$ BCAR-NKT cells (Fig. 4) and compared the gene profiles of these two types of CAR-NKT cell (Fig. 7l). Both populations expressed genes indicative of a mixed T/NK signature (Fig. 7m,n), exhaustion (for example, *BTG1*; Fig. 7o) and cytolytic function (for example, *GZMH*, *GZMK*, *GZMM*, *LAMP1* and *NKG7*; Fig. 7p),

all of which were upregulated after adoptive transfer (Fig. 7m–p). In agreement with the comparison to $PBMC^{15}$ BCAR-NKT cells (Fig. 7a–j), $Allo^{15}$ BCAR-NKT cells demonstrated overall heightened NK gene expression compared with GD2-CAR15 NKT cells (Fig. 7m,n).

Discussion

To generate $Allo^{15}$ CAR-NKT cells compatible with clinical development, we provide a culture system for differentiation and expansion that does not use 3D culture or xenogeneic feeder cells. This technology is versatile, as shown by the successful generation of CAR-NKT cells targeting five different tumor antigens (Fig. 1 and Supplementary Figs. 1, 15 and 16). In addition to CARs, we incorporated other genetic elements into the cells, such as those encoding immune-enhancement molecules (for example, IL-15; Fig. 1d and Supplementary Fig. 1), safety controls (for example, sr39TK; Supplementary Fig. 17)^{18,23} and reporters (for example, Fluc and EGFP; Fig. 3g). None of the evaluated cargo genes appeared to disrupt cell manufacture, yield or quality, suggesting a potential broad application of the technology to target diverse blood cancers, solid tumors and possibly other diseases.

Numerous allogeneic cell therapies have been developed based either on conventional $\alpha\beta$ T cells or NK cells, many equipped with CARs^{5,6,18,69,70}. In the case of the conventional $\alpha\beta$ T cell platform, multiplex gene editing is essential to disrupt the native $\alpha\beta$ TCR, involving the disruption of *TRAC* and/or *TRBC* loci, to prevent GvHD induced by HLA incompatibility⁴. Additionally, HLA-I and/or HLA-II molecules are ablated to mitigate rejection by host T cells⁴, and CD52 is disrupted to render these allogeneic T cells resistant to lymphodepleting drugs such as alemtuzumab⁷¹. The multiplex gene editing increases the risk of aneuploidy and tumorigenicity, and diminishes the cells' efficacy and persistence in vivo^{72,73}. Moreover, ablation of TCR expression compromises cytokine production and long-term persistence compared with CAR-T cells bearing intact TCRs^{74–76}. On the other hand, NK cell-based allogeneic cell products are thought to carry a lower risk of GvHD, simplifying the cell engineering^{5,77}. Nevertheless, their expansion and effectiveness against tumors in vivo might be more restricted compared with conventional $\alpha\beta$ T cells^{5,77}. Various allogeneic cell products have undergone phase I clinical trials targeting B cell malignancies, acute myeloid leukemia, MM, as well as certain solid tumors^{4,5,78,79}, with promising initial results.


A recent phase I clinical trial evaluated an autologous PBMC-derived GD2.CAR-engineered IL-15-enhanced NKT (GD2-CAR15 NKT) cell product in the treatment of relapsed or refractory neuroblastoma^{6,15}. The trial outcomes supported the cells' safety and efficacy^{6,15}. Our finding that $Allo^{15}$ CAR-NKT cells have similar phenotype, functionality and molecular characteristics as their PBMC-derived CAR-NKT cell counterparts underscore the clinical potential of $Allo^{15}$ CAR-NKT cell therapy (Fig. 7). The heightened NK attributes of $Allo^{15}$ CAR-NKT cells compared with PBMC-derived CAR-NKT cells are of interest, particularly for treating solid tumors characterized by loss or heterogeneous expression of the CAR antigen.


Fig. 6 | $Allo^{15}$ CAR-NKT cells exhibit low GvHD risk and high resistance to allojection. **a, b**, Studying the GvH response using an in vitro MLR assay. **a**, Experimental design. **b**, ELISA analyses of IFN γ production on day 4. N, no addition of stimulator PBMCs. **c, d**, Studying the GvHD risk of $Allo^{15}$ BCAR-NKT cells using a human MM xenograft mouse model. Experimental design is shown in Fig. 3a. **c**, Hematoxylin and eosin-stained tissue sections. Tissues were collected on day 60. Scale bars, 100 μ m. **d**, Quantification of **c** ($n = 5$). **e, f**, Studying the T cell-mediated allojection. **e**, Experimental design. **f**, ELISA analyses of IFN γ production on day 4 ($n = 4$). **g, h**, Studying the NK cell-mediated allojection. **g**, Experimental design. **h**, FACS quantification of the indicated live cells ($n = 3$). **i–q**, Studying the molecular mechanisms underlying the allojection resistance of $Allo^{15}$ BCAR-NKT cells. **i**, Illustration depicting the hypoinmunogenicity working model. **j–k**, FACS measurements of surface HLA-I/II and ULBP-1 on $Allo^{15}$ BCAR-NKT cells before infusion (j ; $n = 3$) and post-in vivo antitumor response as described in Fig. 4a (k ; $n = 4$). **l**, Bulk RNA-seq analyses of the cells before infusion ($n = 3$).

m, scRNA-seq analysis of the cells collected from MM xenograft mice on day 28 as described in Fig. 4a. **n**, Methyl-seq analyses of the indicated cells. Vertical lines denote individual CpG sites, wherein red signifies methylation and blue signifies unmethylation. Differentially methylated regions are framed with black boxes. **o**, Western blot analysis of key molecules involved in the IFN γ signaling pathway. **p**, Rescue experiment by overexpressing constitutively active STAT1 (STAT1C) in $Allo^{15}$ BCAR-NKT cells and detecting their surface HLA-I/II and NK ligand expression ($n = 4$). **q**, Schematic illustrating a desensitized IFN γ signaling network within $Allo^{15}$ BCAR-NKT cells. Representative of one (**l–n**) or three (**a–k** and **o–q**) experiments. Data are presented as the mean \pm s.e.m. and were analyzed by two-tailed Student's *t*-test (**d, h, j, k** and **p**) or one-way ANOVA (**b** and **f**). *P* values of violin plots were determined by the Kruskal–Wallis test for the overall comparison, and Dunn's test for post hoc pairwise comparisons between groups (**m**). In the violin plots, box and whisker plots exhibit the minimum, lower quartile, median, upper quartile and maximum expression levels of each type of cell.

Our AlloCAR-NKT cells are free of bystander conventional $\alpha\beta$ T cells (Fig. 1e,g), which reduces the risk of GvHD and enables a streamlined manufacturing process without additional purification steps¹⁶. Because NKT cells do not recognize mismatched HLAs, knockout of

the endogenous TCR is unnecessary, thereby preserving its benefits for CAR-NKT cells⁷⁴. Fifteen cord-blood donors of diverse genetic backgrounds were successfully used to generate CAR-NKT cells (Fig. 1 and Supplementary Fig. 18), avoiding the need for donor selections.





A more than 10⁶-fold expansion was achieved from the input HSPCs to the output mature NKT cells. Assuming a linear extrapolation during the scale-up production, we estimate that from a single cord-blood donor (typically ~1–10 × 10⁶ CD34⁺ HSPCs), roughly 10¹² mature

AlloCAR-NKT cells could be generated, which could potentially be formulated into 1,000–10,000 doses (of ~10⁸–10⁹ cells per dose following current CAR-T cell therapy dosing practices) (Fig. 1a,c and Supplementary Fig. 2c–e)⁸⁰. Note that the high cell yield was only attained when

Fig. 7 | $Allo$ CAR-NKT cells demonstrate similar antitumor capacity to PBMC-derived CAR-NKT cells but with heightened NK attributes. a–j. Comparison between HSPC-derived $Allo^{15}$ BCAR-NKT cells and PBMC-derived IL-15-enhanced BCAR-NKT ($PBMC^{15}$ BCAR-NKT) cells. **a**, Schematics showing the generation of $PBMC^{15}$ BCAR-NKT cells. **b**, FACS plots showing the purity, CAR expression and CD4/CD8 coreceptor expression of $PBMC^{15}$ BCAR-NKT cells. **c**, Yield and CAR expression of $PBMC^{15}$ BCAR-NKT cells ($n = 5$; n indicates different donors). **d**, FACS plots showing the surface expression of NKRs on $Allo^{15}$ BCAR-NKT and $PBMC^{15}$ BCAR-NKT cells. **e**, FACS plots showing the intracellular production of cytokines and cytotoxic molecules in $Allo^{15}$ BCAR-NKT and $PBMC^{15}$ BCAR-NKT cells. **f, g**, Comparing the in vitro antitumor efficacy of $Allo^{15}$ BCAR-NKT and $PBMC^{15}$ BCAR-NKT cells. **f**, Experimental design. **g**, Tumor cell killing data at 24 h ($n = 4$). **h–k**, Comparing the in vivo antitumor efficacy of $Allo^{15}$ BCAR-NKT and $PBMC^{15}$ BCAR-NKT cells in a human MM xenograft NSG mouse model. **h**, Experimental design. **i**, BLI images

showing the tumor loads in experimental mice over time. **j**, Quantification of i ($n = 5$). **k**, Kaplan–Meier survival curves of experimental mice over time ($n = 5$). **l–p**, scRNA-seq comparison between $Allo^{15}$ BCAR-NKT cells and PBMC-derived IL-15-enhanced GD2CAR-engineered NKT (GD2-CAR15 NKT) cells. **l**, Schematics showing the experimental design to compare the gene profiling of $Allo^{15}$ BCAR-NKT and GD2-CAR15 NKT cells using scRNA-seq. GD2-CAR15 NKT cell data are obtained from GSE223071. **m**, Violin plots showing the expression distribution of T and NK cell gene signatures in the indicated CAR-NKT cells. **n–p**, Violin plots showing the expression distribution of genes involved in NK-like differentiation (**n**), exhaustion (**o**) and cytolytic function (**p**), in the indicated CAR-NKT cells. Representative of three (**a–g**), two (**h–k**) and one (**l–p**) experiments. Data are presented as the mean \pm s.e.m. and were analyzed by two-tailed Student's *t*-test (**g**). In the violin plots, box and whisker plots exhibit the minimum, lower quartile, median, upper quartile and maximum expression levels of each type of cell.

culturing TCR-engineered HSPCs but not non-TCR-engineered HSPCs, underscoring a pivotal role of TCR engineering in the reported technology (Supplementary Fig. 2f)^{81,82}.

An additional concern with allogeneic cell products is host cell-mediated allojection. We showed that $Allo$ CAR-NKT cells exhibit an intrinsic, stable hypoimmunogenic phenotype, which is underpinned by epigenetic and signaling network regulation (Fig. 6i–q). Surface HLA-I molecules are only around 20–50% of the levels observed on conventional CAR-T cells (Fig. 6j,k). A recent study of HLA-I knockdown in allogeneic CAR-T cell products⁷⁴ suggests that this level of HLA-I reduction confers resistance against both host T and NK cell-mediated allojection. It is important to note that in vitro MLR assays are models that may not predict allojection in the clinic. Therefore, alternative strategies, such as humanized in vivo models, are needed to better evaluate allojection, although definitive evaluation of allojection can only be accomplished through clinical trials.

Our current $Allo$ CAR-NKT cells comprise mainly CD8 single-positive and/or double-negative cells while missing a CD4 single-positive population that may be of special therapeutic value (Fig. 1e and Supplementary Fig. 1c)^{6,15}; future modification of the culture protocol may generate a CD4 single-positive NKT population. Antitumor potency could be improved using genes encoding other immune-enhancing molecules (for example, IL-7, IL-18 and IL-21) and immunosuppression-resistant factors (for example, immune checkpoint inhibitors such as anti-PD-1 antibody, dominant-negative TGF β receptor)². The safety profile could be further enhanced with suicide switch systems (for example, inducible Cas9 and truncated EGFR)^{83,84}. In vivo performance may also be improved through synergism with other therapeutic modalities (for example, checkpoint blockade therapy and preconditioning regimens)⁷³. Overall, clinical investigation is needed to fully assess the potential of $Allo$ CAR-NKT cells as allogeneic candidates for off-the-shelf cancer therapy.

Online content

Any methods, additional references, Nature Portfolio reporting summaries, source data, extended data, supplementary information, acknowledgements, peer review information; details of author contributions and competing interests; and statements of data and code availability are available at <https://doi.org/10.1038/s41587-024-02226-y>.

References

- June, C. H. & Sadelain, M. Chimeric antigen receptor therapy. *N. Engl. J. Med.* **379**, 64–73 (2018).
- Labanieh, L. & Mackall, C. L. CAR immune cells: design principles, resistance and the next generation. *Nature* **614**, 635–648 (2023).
- Restifo, N. P., Dudley, M. E. & Rosenberg, S. A. Adoptive immunotherapy for cancer: harnessing the T cell response. *Nat. Rev. Immunol.* **12**, 269–281 (2012).
- Benjamin, R. et al. Genome-edited, donor-derived allogeneic anti-CD19 chimeric antigen receptor T cells in paediatric and adult B-cell acute lymphoblastic leukaemia: results of two phase 1 studies. *Lancet* **396**, 1885–1894 (2020).
- Liu, E. et al. Use of CAR-transduced natural killer cells in CD19-positive lymphoid tumors. *N. Engl. J. Med.* **382**, 545–553 (2020).
- Heczey, A. et al. Anti-GD2 CAR-NKT cells in patients with relapsed or refractory neuroblastoma: an interim analysis. *Nat. Med.* **26**, 1686–1690 (2020).
- Klichinsky, M. et al. Human chimeric antigen receptor macrophages for cancer immunotherapy. *Nat. Biotechnol.* **38**, 947–953 (2020).
- Godfrey, D. I., Le Nours, J., Andrews, D. M., Uldrich, A. P. & Rossjohn, J. Unconventional T cell targets for cancer immunotherapy. *Immunity* **48**, 453–473 (2018).
- Chandra, S. & Kronenberg, M. Activation and function of iNKT and MAIT cells. *Adv. Immunol.* **127**, 145–201 (2015).
- Chaidos, A. et al. Graft invariant natural killer T-cell dose predicts risk of acute graft-versus-host disease in allogeneic hematopoietic stem cell transplantation. *Blood* **119**, 5030–5036 (2012).
- Exley, M. A. et al. Adoptive transfer of invariant NKT cells as immunotherapy for advanced melanoma: a phase I clinical trial. *Clin. Cancer Res.* **23**, 3510–3519 (2017).
- Fang, Y. et al. Graft-versus-host disease modulation by innate T cells. *Int. J. Mol. Sci.* **24**, 4084 (2023).
- Shlomchik, W. D. Graft-versus-host disease. *Nat. Rev. Immunol.* **7**, 340–352 (2007).
- Brennan, P. J., Brigl, M. & Brenner, M. B. Invariant natural killer T cells: an innate activation scheme linked to diverse effector functions. *Nat. Rev. Immunol.* **13**, 101–117 (2013).
- Heczey, A. et al. Anti-GD2 CAR-NKT cells in relapsed or refractory neuroblastoma: updated phase 1 trial interim results. *Nat. Med.* <https://doi.org/10.1038/s41591-023-02363-y> (2023).
- Li, Y.-R. et al. Advancing cell-based cancer immunotherapy through stem cell engineering. *Cell Stem Cell* <https://doi.org/10.1016/j.stem.2023.02.009> (2023).
- Gschwend, E., De Oliveira, S. & Kohn, D. B. Hematopoietic stem cells for cancer immunotherapy. *Immunol. Rev.* **257**, 237–249 (2014).
- Li, Y.-R. et al. Development of allogeneic HSC-engineered iNKT cells for off-the-shelf cancer immunotherapy. *Cell Rep. Med.* **2**, 100449 (2021).
- Smith, D. J. et al. Genetic engineering of hematopoietic stem cells to generate invariant natural killer T cells. *Proc. Natl Acad. Sci. USA* **112**, 1523–1528 (2015).
- Pilipow, K. et al. IL15 and T-cell stemness in T-cell-based cancer immunotherapy. *Cancer Res.* **75**, 5187–5193 (2015).

21. Zhou, Y. et al. Interleukin 15 in Cell-Based Cancer Immunotherapy. *Int. J. Mol. Sci.* **23**, 7311 (2022).
22. Gordy, L. E. et al. IL-15 regulates homeostasis and terminal maturation of NKT cells. *J. Immunol.* **187**, 6335–6345 (2011).
23. Zhu, Y. et al. Development of hematopoietic stem cell-engineered invariant natural killer T cell therapy for cancer. *Cell Stem Cell* **25**, 542–557.e9 (2019).
24. Liu, Y. et al. Durable remission achieved from BCMA-directed CAR-T therapy against relapsed or refractory multiple myeloma. *Blood* **132**, 956 (2018).
25. Cheng, Z. et al. In vivo expansion and antitumor activity of coinfused CD28- and 4-1BB-engineered CAR-T cells in patients with B cell leukemia. *Mol. Ther.* **26**, 976–985 (2018).
26. Heczey, A. et al. Invariant NKT cells with chimeric antigen receptor provide a novel platform for safe and effective cancer immunotherapy. *Blood* **124**, 2824–2833 (2014).
27. Jiang, Z. et al. Anti-GPC3-CAR T cells suppress the growth of tumor cells in patient-derived xenografts of hepatocellular carcinoma. *Front. Immunol.* **7**, 690 (2016).
28. O'Rourke, D. M. et al. A single dose of peripherally infused EGFRvIII-directed CAR T cells mediates antigen loss and induces adaptive resistance in patients with recurrent glioblastoma. *Sci. Transl. Med.* **9**, eaaa0984 (2017).
29. Hurton, L. V. et al. Tethered IL-15 augments antitumor activity and promotes a stem-cell memory subset in tumor-specific T cells. *Proc. Natl Acad. Sci. USA* **113**, E7788–E7797 (2016).
30. Xu, X. et al. NKT cells coexpressing a GD2-specific chimeric antigen receptor and IL15 show enhanced in vivo persistence and antitumor activity against neuroblastoma. *Clin. Cancer Res.* **25**, 7126–7138 (2019).
31. Li, Y.-R., Zhou, Y., Kramer, A. & Yang, L. Engineering stem cells for cancer immunotherapy. *Trends Cancer* <https://doi.org/10.1016/j.trecan.2021.08.004> (2021).
32. Godfrey, D. I. & Berzins, S. P. Control points in NKT-cell development. *Nat. Rev. Immunol.* **7**, 505–518 (2007).
33. Giannoni, F. et al. Allelic exclusion and peripheral reconstitution by TCR transgenic T cells arising from transduced human hematopoietic stem/progenitor cells. *Mol. Ther.* **21**, 1044–1054 (2013).
34. Lam, P. Y., Nissen, M. D. & Mattarollo, S. R. Invariant natural killer T cells in immune regulation of blood cancers: harnessing their potential in immunotherapies. *Front. Immunol.* **8**, 1355 (2017).
35. Krijgsman, D., Hokland, M. & Kuppen, P. J. K. The role of natural killer T cells in cancer—a phenotypical and functional approach. *Front. Immunol.* **9**, 367 (2018).
36. Li, Y.-R. et al. Off-the-shelf third-party HSC-engineered iNKT cells for ameliorating GvHD while preserving GvL effect in the treatment of blood cancers. *iScience* **25**, 104859 (2022).
37. Sommer, C. et al. Preclinical evaluation of allogeneic CAR T cells targeting BCMA for the treatment of multiple myeloma. *Mol. Ther.* **27**, 1126–1138 (2019).
38. Miller, B. C. et al. Subsets of exhausted CD8(+) T cells differentially mediate tumor control and respond to checkpoint blockade. *Nat. Immunol.* **20**, 326–336 (2019).
39. Zheng, L. et al. Pan-cancer single-cell landscape of tumor-infiltrating T cells. *Science* **374**, abe6474 (2021).
40. Li, H. et al. Dysfunctional CD8 T cells form a proliferative, dynamically regulated compartment within human melanoma. *Cell* **176**, 775–789.e18 (2019).
41. Milner, J. J. et al. Heterogenous populations of tissue-resident CD8(+) T cells are generated in response to infection and malignancy. *Immunity* **52**, 808–824.e7 (2020).
42. Beltra, J.-C. et al. Developmental relationships of four exhausted CD8(+) T cell subsets reveals underlying transcriptional and epigenetic landscape control mechanisms. *Immunity* **52**, 825–841.e8 (2020).
43. Zheng, C. et al. Landscape of infiltrating T cells in liver cancer revealed by single-cell sequencing. *Cell* **169**, 1342–1356.e16 (2017).
44. Wang, X., Peticone, C., Kotsopoulou, E., Göttgens, B. & Calero-Nieto, F. J. Single-cell transcriptome analysis of CAR T-cell products reveals subpopulations, stimulation, and exhaustion signatures. *Oncoimmunology* **10**, 1866287 (2021).
45. Bai, Z. et al. Single-cell antigen-specific landscape of CAR T infusion product identifies determinants of CD19-positive relapse in patients with ALL. *Sci. Adv.* **8**, eabj2820 (2022).
46. Good, C. R. et al. An NK-like CAR T cell transition in CAR T cell dysfunction. *Cell* **184**, 6081–6100.e26 (2021).
47. Haradhvala, N. J. et al. Distinct cellular dynamics associated with response to CAR-T therapy for refractory B cell lymphoma. *Nat. Med.* **28**, 1848–1859 (2022).
48. Guo, X. et al. Global characterization of T cells in non-small-cell lung cancer by single-cell sequencing. *Nat. Med.* **24**, 978–985 (2018).
49. Liu, D. et al. IL-15 protects NKT cells from inhibition by tumor-associated macrophages and enhances antimetastatic activity. *J. Clin. Invest.* **122**, 2221–2233 (2012).
50. Rafiq, S., Hackett, C. S. & Brentjens, R. J. Engineering strategies to overcome the current roadblocks in CAR T cell therapy. *Nat. Rev. Clin. Oncol.* **17**, 147–167 (2020).
51. Majzner, R. G. & Mackall, C. L. Tumor antigen escape from CAR T-cell therapy. *Cancer Discov.* **8**, 1219–1226 (2018).
52. Li, Y.-R. et al. Targeting immunosuppressive tumor-associated macrophages using innate T cells for enhanced antitumor reactivity. *Cancers* **14**, 2749 (2022).
53. Li, Y.-R. et al. An ex vivo 3D tumor microenvironment-mimicry culture to study TAM modulation of cancer immunotherapy. *Cells* **11**, 1583 (2022).
54. Li, Y.-R., Wilson, M. & Yang, L. Target tumor microenvironment by innate T cells. *Front. Immunol.* **13**, 999549 (2022).
55. Matsuda, J. L., Mallevaey, T., Scott-Browne, J. & Gapin, L. CD1d-restricted iNKT cells, the ‘Swiss-Army knife’ of the immune system. *Curr. Opin. Immunol.* **20**, 358–368 (2008).
56. Kawano, Y. et al. Targeting the bone marrow microenvironment in multiple myeloma. *Immunol. Rev.* **263**, 160–172 (2015).
57. Neumeister, P., Schulz, E., Pansy, K., Szmyra, M. & Deutsch, A. J. A. Targeting the microenvironment for treating multiple myeloma. *Int. J. Mol. Sci.* **23**, 7627 (2022).
58. Sun, J. et al. Tumor-associated macrophages in multiple myeloma: advances in biology and therapy. *J. Immunother. Cancer* **10**, e003975 (2022).
59. Giavridis, T. et al. CAR T cell-induced cytokine release syndrome is mediated by macrophages and abated by IL-1 blockade letter. *Nat. Med.* **24**, 731–738 (2018).
60. Norelli, M. et al. Monocyte-derived IL-1 and IL-6 are differentially required for cytokine-release syndrome and neurotoxicity due to CAR T cells. *Nat. Med.* **24**, 739–748 (2018).
61. Suárez-Alvarez, B. et al. Epigenetic mechanisms regulate MHC and antigen processing molecules in human embryonic and induced pluripotent stem cells. *PLoS ONE* **5**, e10192 (2010).
62. Majumder, P. & Boss, J. M. DNA methylation dysregulates and silences the HLA-DQ locus by altering chromatin architecture. *Genes Immun.* **12**, 291–299 (2011).
63. Liu, X. et al. Activation of the JAK/STAT-1 signaling pathway by IFN-gamma can down-regulate functional expression of the MHC class I-related neonatal Fc receptor for IgG. *J. Immunol.* **181**, 449–463 (2008).
64. Hu, X., Li, J., Fu, M., Zhao, X. & Wang, W. The JAK/STAT signaling pathway: from bench to clinic. *Signal Transduct. Target. Ther.* **6**, 402 (2021).

65. Jongsma, M. L. M., Guarda, G. & Spaapen, R. M. The regulatory network behind MHC class I expression. *Mol. Immunol.* **113**, 16–21 (2019).
66. Masternak, K. et al. CIITA is a transcriptional coactivator that is recruited to MHC class II promoters by multiple synergistic interactions with an enhanceosome complex. *Genes Dev.* **14**, 1156–1166 (2000).
67. Tan, N. Y. & Khachigian, L. M. Sp1 phosphorylation and its regulation of gene transcription. *Mol. Cell. Biol.* **29**, 2483–2488 (2009).
68. Liddle, F. J., Alvarez, J. V., Poli, V. & Frank, D. A. Tyrosine phosphorylation is required for functional activation of disulfide-containing constitutively active STAT mutants. *Biochemistry* **45**, 5599–5605 (2006).
69. van der Stegen, S. J. C. et al. Generation of T-cell-receptor-negative CD8 α β -positive CAR T cells from T-cell-derived induced pluripotent stem cells. *Nat. Biomed. Eng.* **6**, 1284–1297 (2022).
70. Wang, Z. et al. 3D-organoid culture supports differentiation of human CAR+ iPSCs into highly functional CAR T cells. *Cell Stem Cell* **29**, 515–527.e8 (2022).
71. Benjamin, R. et al. Preliminary data on safety, cellular kinetics and anti-leukemic activity of UCART19, an allogeneic anti-CD19 CAR T-cell product, in a pool of adult and pediatric patients with high-risk CD19+ relapsed/refractory B-cell acute lymphoblastic leukemia. *Blood* **132**, 896–896 (2018).
72. Nahmad, A. D. et al. Frequent aneuploidy in primary human T cells after CRISPR-Cas9 cleavage. *Nat. Biotechnol.* **40**, 1807–1813 (2022).
73. Depil, S., Duchateau, P., Grupp, S. A., Mufti, G. & Poirot, L. Off-the-shelf' allogeneic CAR T cells: development and challenges. *Nat. Rev. Drug Discov.* **19**, 185–199 (2020).
74. Xia, H. et al. Presidential symposium and presentation of top abstracts. *Mol. Ther.* **31**, 1–794 (2023).
75. Stenger, D. et al. Endogenous TCR promotes in vivo persistence of CD19-CAR-T cells compared to a CRISPR/Cas9-mediated TCR knockout CAR. *Blood* **136**, 1407–1418 (2020).
76. Wang, Z. et al. Phase I study of CAR-T cells with PD-1 and TCR disruption in mesothelin-positive solid tumors. *Cell. Mol. Immunol.* **18**, 2188–2198 (2021).
77. Cichocki, F. et al. iPSC-derived NK cells maintain high cytotoxicity and enhance in vivo tumor control in concert with T cells and anti-PD-1 therapy. *Sci. Transl. Med.* **12**, eaaz5618 (2020).
78. Benjamin, R. et al. UCART19, a first-in-class allogeneic anti-CD19 chimeric antigen receptor T-cell therapy for adults with relapsed or refractory B-cell acute lymphoblastic leukaemia (CALM): a phase 1, dose-escalation trial. *Lancet Haematol.* **9**, e833–e843 (2022).
79. Mailankody, S. et al. Allogeneic BCMA-targeting CAR T cells in relapsed/refractory multiple myeloma: phase 1 UNIVERSAL trial interim results. *Nat. Med.* **29**, 422–429 (2023).
80. Labanieh, L., Majzner, R. G. & Mackall, C. L. Programming CAR-T cells to kill cancer. *Nat. Biomed. Eng.* **2**, 377–391 (2018).
81. Ellis, G. I., Sheppard, N. C. & Riley, J. L. Genetic engineering of T cells for immunotherapy. *Nat. Rev. Genet.* **22**, 427–447 (2021).
82. Oppermans, N., Kueberuwa, G., Hawkins, R. E. & Bridgeman, J. S. Transgenic T-cell receptor immunotherapy for cancer: building on clinical success. *Ther. Adv. Vaccines Immunother.* **8**, 2515135520933509 (2020).
83. Di Stasi, A. et al. Inducible apoptosis as a safety switch for adoptive cell therapy. *N. Engl. J. Med.* **365**, 1673–1683 (2011).
84. Straathof, K. C. et al. An inducible caspase 9 safety switch for T-cell therapy. *Blood* **105**, 4247–4254 (2005).

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>.

© The Author(s) 2024

¹Department of Microbiology, Immunology and Molecular Genetics, University of California, Los Angeles, Los Angeles, CA, USA. ²Department of Medicine, Division of Cardiology, University of California, Los Angeles, Los Angeles, CA, USA. ³Mork Family Department of Chemical Engineering and Materials Science, University of Southern California, Los Angeles, CA, USA. ⁴Bioinformatics Interdepartmental Program, University of California, Los Angeles, Los Angeles, CA, USA. ⁵Department of Pathology and Laboratory Medicine, University of California, Los Angeles, Los Angeles, CA, USA. ⁶Department of Cell and Developmental Biology, Feinberg School of Medicine, Northwestern University, Chicago, IL, USA. ⁷Human Genetics, University of California, Los Angeles, Los Angeles, CA, USA. ⁸Department of Biostatistics, Fielding School of Public Health, University of California, Los Angeles, Los Angeles, CA, USA. ⁹Eli and Edythe Broad Centre of Regenerative Medicine and Stem Cell Research, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁰Jonsson Comprehensive Cancer Centre, University of California, Los Angeles, Los Angeles, CA, USA. ¹¹Department of Medicine, Division of Hematology/Oncology, University of California, Los Angeles, Los Angeles, CA, USA. ¹²Molecular Biology Institute, University of California, Los Angeles, Los Angeles, CA, USA. ¹³Department of Pediatrics, Division of Hematology/Oncology, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁴Institute for Quantitative and Computational Biosciences—The Collaboratory, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁵Department of Molecular, Cell, and Developmental Biology, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁶Department of Hematology and Oncology, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁷Department of Internal Medicine, University of California, Los Angeles, Los Angeles, CA, USA. ¹⁸These authors contributed equally: Yan-Ruide Li, Yang Zhou, Jiaji Yu. ✉ e-mail: liliyang@ucla.edu

Methods

Mice

NOD.Cg-Prkdc^{SCID}Il2rg^{tm1Wjl}/SzJ (NOD/SCID/IL-2Rγ^{-/-}, NSG) mice were maintained in the animal facilities of University of California, Los Angeles (UCLA). Here 6–10 week-old female mice were used for all experiments unless otherwise indicated. All animal experiments were approved by the Institutional Animal Care and Use Committee of UCLA. All mice were bred and maintained under specific pathogen-free conditions, and all experiments were conducted in accordance with the animal care and use regulations of the Division of Laboratory Animal Medicine at the UCLA.

Media and reagents

The X-VIVO 15 Serum-Free Hematopoietic Cell Medium was purchased from Lonza. The StemSpan T Cell Generation Kit, comprising the StemSpan SFEM II Medium, the StemSpan Lymphoid Progenitor Expansion Supplement, the StemSpan Lymphoid Progenitor Maturation Supplement, the StemSpan Lymphoid Progenitor Differentiation Coating Material and the ImmunoCult Human CD3/CD28/CD2 T Cell Activator, was purchased from StemCell Technologies. The CTS OpTmizer T Cell Expansion SFM (no phenol red, bottle format), the Roswell Park Memorial Institute (RPMI) 1640 cell culture medium and the DMEM cell culture medium were purchased from Thermo Fisher Scientific. The CryoStor Cell Cryopreservation Media CS10 was purchased from MilliporeSigma.

α-Galactosylceramide (αGC, KRN7000) was purchased from Avanti Polar Lipids. Recombinant human IL-2, IL-3, IL-7, IL-15, IL-21, IFNγ, Flt3 ligand (Flt3L), stem cell factor (SCF) and thrombopoietin (TPO) were purchased from Peprotech. Ganciclovir, fetal bovine serum (FBS) and beta-mercaptoethanol (β-ME) were purchased from Sigma. Penicillin-streptomycin-glutamine (PSG), MEM non-essential amino acids, HEPES buffer solution and sodium pyruvate were purchased from Gibco. Normocin was purchased from Invivogen.

The homemade C10 medium was made of RPMI 1640 cell culture medium, supplemented with FBS (10% vol/vol), PSG (1% vol/vol), MEM non-essential amino acids (1% vol/vol), HEPES (10 mM), sodium pyruvate (1 mM), β-ME (50 mM) and normocin (100 mg ml⁻¹). The homemade D10 medium was made of DMEM supplemented with FBS (10% vol/vol), PSG (1% vol/vol) and normocin (100 mg ml⁻¹). The homemade R10 medium was made of RPMI supplemented with FBS (10% vol/vol), PSG (1% vol/vol) and normocin (100 mg ml⁻¹).

Lentiviral vectors

Lentiviral vectors used in this study were all constructed from a parental lentivector pMNDW^{18,23,33}. The 2A sequences derived from foot-and-mouth disease virus (F2A), porcine teschovirus-1 (P2A) and theosea asigna virus (T2A) were used to link the inserted genes to achieve coexpression.

The Lenti/iNKT vector was constructed by inserting into the pMNDW parental vector a synthetic bicistronic gene encoding human iNKT TCRα-F2A-TCRβ (refs. 18,23). The Lenti/iNKT-sr39TK vector was constructed by inserting a synthetic tricistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-sr39TK. The Lenti/iNKT-BCAR vector was constructed by inserting a synthetic tricistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-BCAR (BCAR indicates a BCMA-targeting CAR)²⁴. The Lenti/iNKT-BCAR-IL15 vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-BCAR-T2A-IL15 (IL15 indicates the secreted form of human IL-15). The Lenti/iNKT-CAR19 vector was constructed by inserting a synthetic tricistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-CAR19 (CAR19 indicates a CD19-targeting CAR)²⁵. The Lenti/iNKT-CAR19-IL15 vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-CAR19-T2A-IL15. The Lenti/iNKT-GD2.CAR vector was constructed by inserting a synthetic tricistronic gene encoding

human iNKT TCRα-F2A-TCRβ-P2A-GD2.CAR (GD2.CAR indicates a GD2-targeting CAR)²⁶. The Lenti/iNKT-GD2.CAR-IL15 vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-GD2.CAR-T2A-IL15. The Lenti/iNKT-GPC3.CAR vector was constructed by inserting a synthetic tricistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-GPC3.CAR (GPC3.CAR indicates a glypican-3-targeting CAR)²⁷. The Lenti/iNKT-GPC3.CAR-IL15 vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-GPC3.CAR-T2A-IL15. The Lenti/iNKT-ECAR vector was constructed by inserting a synthetic tricistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-ECAR (ECAR indicates an epidermal growth factor receptor variant III-targeting CAR)²⁸. The Lenti/iNKT-ECAR-IL15 vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-ECAR-T2A-IL15. The Lenti/iNKT-IL15-sr39TK vector was constructed by inserting a synthetic tetracistronic gene encoding human iNKT TCRα-F2A-TCRβ-P2A-IL15-T2A-sr39TK (sr39TK indicates an sr39TK suicide and positron emission tomography imaging reporter gene).

The Lenti/FlucGFP vector was constructed by inserting into the pMNDW a synthetic bicistronic gene encoding Fluc-P2A-EGFP (ref. 23). The Lenti/CD1d vector was constructed by inserting a synthetic gene encoding human CD1d (ref. 23). The Lenti/EGFRvIII vector was constructed by inserting a synthetic gene encoding human EGFRvIII. The Lenti/STAT1C-EGFP vector was constructed by inserting a synthetic bicistronic gene encoding STAT1C-P2A-EGFP (STAT1C indicates a human constitutively active STAT1, ref. 85). The Lenti/BCAR vector was constructed by inserting a synthetic gene encoding BCAR. The Lenti/BCAR-IL15 vector was constructed by inserting a synthetic bicistronic gene encoding BCAR-F2A-IL15. The Lenti/CAR19 vector was constructed by inserting a synthetic gene encoding CAR19. The Lenti/GD2.CAR vector was constructed by inserting a synthetic gene encoding GD2.CAR. The Lenti/GPC3.CAR vector was constructed by inserting a synthetic gene encoding GPC3.CAR. The Lenti/ECAR vector was constructed by inserting a synthetic gene encoding ECAR.

The synthetic gene fragments were obtained from GenScript and IDT. Lentiviruses were produced using human embryonic kidney (HEK) 293T cells (American Type Culture Collection (ATCC)), following a standard transfection protocol using the Trans-IT-Lenti Transfection Reagent (Mirus Bio) and a centrifugation concentration protocol using the Amicon Ultra Centrifugal Filter Units, according to the manufacturers' instructions (MilliporeSigma).

Cell lines

Human MM cell line MM.1S, chronic myelogenous leukemia cell line K562, Burkitt's lymphoma cell line Raji, acute lymphoblastic leukemia cell line NALM-6, melanoma cell line A375, glioblastoma cell line T98G and U87MG, hepatocellular carcinoma cell line HEP3B and HEK293T were purchased from the ATCC.

To make stable tumor cell lines overexpressing human CD1d or EGFRvIII, and/or FG dual reporters, the parental tumor cell lines were transduced with lentiviral vectors encoding the intended gene(s). Then 72 h after lentivector transduction, cells were subjected to flow cytometry sorting to isolate gene-engineered cells for making stable cell lines. Nine stable tumor cell lines were generated for this study, including MM-FG, MM-FG-CD1d, K562-FG, Raji-FG, Raji-FG-CD1d, NALM-6-FG, T98G-FG, U87MG-EGFRvIII-FG and Hep3B-FG cell lines. The ^{KO}MM-FG cell line was generated by knocking out the BCMA gene from the parental MM-FG cell line using CRISPR-Cas9. The single-guide RNA (sgRNA) targeting the BCMA gene (UAUUAAGCUCAGUCCCAAAC) was purchased from Synthego, and was introduced into MM-FG cells via electroporation using an Amaxa 4D Nucleofection X Unit (Lonza), according to the manufacturer's instructions.

The artificial antigen-presenting cell (aAPC) was generated by engineering the K562 cell line to overexpress human CD83/CD86/4-1BBL costimulatory receptors⁸⁶. The aAPC-BCMA, aAPC-CD19,

aAPC-GD2, aAPC-GPC3 and aAPC-EGFRvIII cell lines were generated by further engineering the parental aAPC line to overexpress human BCMA, CD19, GD2, GPC3 and EGFRvIII, respectively.

Human CD34⁺ HSPCs, PBMCs and primary BM samples from patients with MM

Purified cord-blood-derived human CD34⁺ cells were purchased from the HemaCare. Healthy donor PBMCs were obtained from the UCLA/CFAR Virology Core Laboratory without identification information under federal and state regulations. Primary samples from patients with MM were collected at the Ronald Reagan UCLA Medical Center from consented patients through an IRB-approved protocol (IRB no. 21-001444) and processed.

Antibodies and flow cytometry

Fluorochrome-conjugated antibodies specific for human CD45 (Clone HI30, cat. no. 304026, 1:500 dilution), TCR $\alpha\beta$ (Clone I26, cat. no. 306716, 1:25 dilution), CD3 (Clone HIT3a, cat. no. 300329, 1:500 dilution), CD4 (Clone OKT4, cat. no. 317414, 1:400 dilution), CD8 (Clone SK1, cat. no. 344714, 1:500 dilution), CD45RO (Clone UCHL1, cat. no. 304216, 1:200 dilution), CD45RA (Clone HI100, cat. no. 304105, 1:5,000 dilution), CD161 (Clone HP-3G10, cat. no. 339928, 1:50 dilution), CD69 (Clone FN50, cat. no. 310909, 1:50 dilution), CD56 (Clone HCD56, cat. no. 362545, 1:10 dilution), CD1d (Clone 51.1, cat. no. 350308, 1:50 dilution), BCMA (19F2, cat. no. 357503, 1:50 dilution), CD14 (Clone HCD14, cat. no. 325608, 1:100 dilution), CD19 (Clone HIB19, cat. no. 363005, 1:100 dilution), CD11b (Clone ICRF44, cat. no. 301330, 1:500 dilution), CD15 (W6D3, cat. no. 323021, 1:500 dilution), CD112 (Clone TX31, cat. no. 337409, 1:200 dilution), CD155 (Clone SKII.4, cat. no. 337613, 1:200 dilution), MICA/MICB (Clone 6D4, cat. no. 320908, 1:50 dilution), Ganglioside GD2 (14G2a, cat. no. 357323, 1:50 dilution), NKG2D (Clone 1D11, cat. no. 320812, 1:50 dilution), DNAM-1 (Clone 11A8, cat. no. 338312, 1:50 dilution), CD158 (KIR2DL1/S1/S3/S5) (Clone HP-MA4, cat. no. 339510, 1:50 dilution), Nkp30 (Clone P30-15, cat. no. 325207, 1:50 dilution), Nkp44 (Clone P44-8, cat. no. 325107, 1:50 dilution), CD16 (Clone 3G8, cat. no. 302011, 1:50 dilution), NKG2A (Clone S19004C, cat. no. 375103, 1:50 dilution), NKG2C (Clone S19005E, cat. no. 375003, 1:50 dilution), CD62L (Clone P44-8, cat. no. 304813, 1:50 dilution), CD134 (Clone Ber-ACT35, cat. no. 350008, 1:50 dilution), LEF1 (Clone W17021C, cat. no. 621051, 1:50 dilution), CD8A (Clone C8/144B, cat. no. 372092, 1:500 dilution), CD8B (Clone QA20A40, cat. no. 387305, 1:500 dilution), IFN γ (Clone B27, cat. no. 506518, 1:50 dilution), Granzyme B (Clone QA16A02, cat. no. 372204, 1:4,000 dilution), Perforin (Clone dG9, cat. no. 308126, 1:50 dilution), TNF α (Clone Mab11, cat. no. 502912, 1:4,000 dilution), IL-2 (Clone MQ1-17H12, cat. no. 500341, 1:200 dilution), β 2-microglobulin (B2M) (Clone 2M2, cat. no. 316312, 1:5,000 dilution), HLA-DR (Clone L243, cat. no. 307618, 1:250 dilution), and HLA-DR, double-positive, DQ (Clone Tü 39, cat. no. 361707, 1:250 dilution) were purchased from BioLegend. Fluorochrome-conjugated antibodies specific for mouse CD45 (Clone S18009F, cat. no. 157607, 1:5,000 dilution), GR1 (Clone RB6-8C5, cat. no. 108411, 1:1,000 dilution), CD1d (Clone IB1, cat. no. 123521, 1:100 dilution), CD11b (Clone MI/70, cat. no. 101205, 1:5,000 dilution), Sca1 (Clone D7, cat. no. 108111, 1:100 dilution), Flt3 (Clone A2F10, cat. no. 135305, 1:100 dilution), SLAM (Clone TC15-12F12.2, cat. no. 115913, 1:50 dilution), c-Kit (Clone 2B8, cat. no. 105825, 1:50 dilution) were purchased from BioLegend. In our study, note the use of antibodies with identical clones but differing conjugated fluorochromes, with one typical antibody listed herein. Fluorochrome-conjugated antibody specific for human Glypican-3 (GPC3; Clone 024, cat. no. ab275695, 1:200 dilution) was purchased from Abcam. Fluorochrome-conjugated antibodies specific for human CD34 (Clone 581, cat. no. 555822, 1:100 dilution) and human iNKT TCR V α 24-J β 18 (Clone 6B11, cat. no. 552825, 1:10 dilution) were purchased from BD Biosciences. Fluorochrome-conjugated antibody specific for human iNKT TCR V β 11 (Clone C21, cat. no. A66905, 1:50 dilution)

was purchased from Beckman-Coulter. Fluorochrome-conjugated antibodies specific for human ULBP-1 (Clone 170818, cat. no. FAB1380P, 1:50 dilution) and ULBP-2,5,6 (Clone I65903, cat. no. FAB1298A5, 1:50 dilution) were purchased from R&D Systems. A goat antimouse IgG F(ab')₂ secondary antibody (cat. no. 31803, 1:50 dilution) was purchased from Thermo Fisher. Fixable Viability Dye eFluor506 (e506, cat. no. 65-0866-14, 1:500 dilution) was purchased from Affymetrix eBioscience; mouse Fc Block (antimouse CD16/32, cat. no. 553142, 1:50 dilution) was purchased from BD Biosciences and human Fc Receptor Blocking Solution (TrueStain FcX, cat. no. 422302, 1:100 dilution) was purchased from BioLegend.

All flow cytometry staining was performed following standard protocols, as well as specific instructions provided by the manufacturer of a particular antibody. Stained cells were analyzed using a MACSQuant Analyzer 10 flow cytometer (Miltenyi Biotech), following the manufacturer's instructions. FlowJo software v.9 (BD Biosciences) was used for data analysis.

Enzyme-linked immunosorbent cytokine assays (ELISAs)

Supernatants from cell culture assays were collected and assayed to quantify human IFN γ , TNF α , IL-2, and IL-4. The capture and biotinylated pairs for detecting cytokines were purchased from BD Biosciences. The streptavidin-HRP conjugate was purchased from Invitrogen. Human cytokine standards were purchased from eBioscience. Tetramethylbenzidine (TMB) substrate was purchased from KPL. Human IL-15 was quantified using a Human IL-15 Quantikine ELISA Kit (R&D Systems), following the manufacturer's instructions. Human IL-17a was quantified using a Human IL-17A ELISA MAX Deluxe Kit (BioLegend), following the manufacturer's instructions. Mouse IL-6 was quantified with paired purified anti-mouse IL-6 antibody and biotin anti-mouse IL-6 antibody (BioLegend). Mouse SAA-3 was quantified using a Mouse SAA-3 ELISA Kit (MilliporeSigma), per the manufacturer's instructions. The samples were analyzed for absorbance at 450 nm using an Infinite M1000 microplate reader (Tecan).

Generation of allogeneic HSPC-engineered NKT (Allo¹⁵NKT) cells and their CAR/IL15-armed derivatives (denoted as Allo¹⁵CAR-NKT cells)

Allo¹⁵NKT and Allo¹⁵CAR-NKT cells were generated by differentiating gene-engineered cord-blood CD34⁺ HSPCs in a five-stage Ex Vivo HSPC-Derived NKT Cell Culture. Allo¹⁵NKT cells were differentiated from HSPCs engineered to overexpress a transgenic human iNKT TCR, while Allo¹⁵CAR-NKT cells were differentiated from HSPCs engineered to overexpress a human transgenic iNKT TCR, together with a selected CAR and/or the secreted form of human IL-15.

At stage 0, 1×10^4 frozen-thawed human CD34⁺ HSPCs were revived and cultured in 300 μ l of X-VIVO 15 Serum-Free Hematopoietic Stem Cell Medium supplemented with 50 ng ml⁻¹ Flt3L, 50 ng ml⁻¹ stem cell factor, 50 ng ml⁻¹ thrombopoietin and 20 ng ml⁻¹ IL-3 in non-tissue culture treated 24-well plate for 24 h, then transduced with Lenti/iNKT-(CAR)-(IL15) viruses for another 24 h following an established protocol^{18,23}. Briefly, concentrated lentivirus supernatant was mixed with 1:100 Poloxamer Synperonic F108 and 1:1,000 Prostaglandin E2 (PGE2; CAYMAN), and then gently added to the HSPCs culture.

At stage 1, gene-engineered HSPCs collected from stage 0 were cultured in the StemSpan SFEM II Medium supplemented with StemSpan Lymphoid Progenitor Expansion Supplement for 2 weeks. CELLSTAR 24-well Cell Culture Non-treated Multiwell Plates (VWR) were used. The plates were coated with 500 μ l per well StemSpan Lymphoid Differentiation Coating Material for 2 h at room temperature or alternatively, overnight at 4 °C. Transduced CD34⁺ HSPCs were suspended at 2×10^4 cells per ml and 500 μ l of cell suspension was added into each precoated well. Twice per week, half of the medium from each well was removed and replaced with fresh medium.

At stage 2, cells collected from the stage 1 were cultured in the StemSpan SFEM II Medium supplemented with StemSpan Lymphoid Progenitor Maturation Supplement for 1 week. Non-Treated Falcon Polystyrene six-well Microplates (Thermo Fisher Scientific) were coated with 1 ml per well of StemSpan Lymphoid Differentiation Coating Material. The stage 1 cells were collected and resuspended at 1×10^5 cells per ml; 2 ml of cell suspension was added into each precoated well. Cells were passaged 1–2 times per week to maintain a cell density at $0.5\text{--}1 \times 10^6$ cells per ml; fresh medium was added at every passage.

At stage 3, cells collected from the stage 2 were cultured in the StemSpan SFEM II Medium supplemented with StemSpan Lymphoid Progenitor Maturation Supplement, CD3/CD28/CD2 T Cell Activator and 20 ng ml⁻¹ human recombinant IL-15 for 1 week. Cells were resuspended at 5×10^5 cells per ml; 2 ml cell suspension was added into Non-Treated Falcon Polystyrene six-well Microplates (Thermo Fisher Scientific) precoated with 1 ml per well of StemSpan Lymphoid Differentiation Coating Material. Cells were passaged 2–3 times per week to maintain a cell density at $0.5\text{--}1 \times 10^6$ cells per ml; fresh medium was added at every passage.

At stage 4, cells collected from the stage 3, now mature Allo¹⁵(CAR)-NKT cells or their derivatives, were expanded using various expansion approaches: (1) an α CD3/ α CD28 expansion approach, (2) an α GC-PBMC expansion approach or (3) an aAPC expansion approach. The expansion stage lasted for 2 weeks. At stage 4, cells could be cultured in 150 mm cell culture dishes (Thermo Fisher Scientific) or G-Rex 6M Well Plates (Wilson Wolf). The expansion can happen in a feeder-free, serum-free CTS OpTmizer T Cell Expansion SFM (Thermo Fisher Scientific) or a homemade C10 medium. The resulting Allo¹⁵(CAR)-NKT or derivative cell products were aliquoted and cryopreserved in CryoStor Cell Cryopreservation Media CS10 using a Thermo Scientific CryoMed Controlled-Rate Freezer 7450 (Thermo Scientific) for future use, following the manufacturer's instructions.

- (1) The α CD3/ α CD28 antibody expansion approach used 150 mm cell culture dishes (Thermo Fisher Scientific), which were coated with 1 μ g ml⁻¹ (500 μ l per well) of Ultra-LEAF Purified Anti-Human CD3 Antibody (Clone OKT3; BioLegend) for 2 h at room temperature or, alternatively, overnight at 4 °C. Mature Allo¹⁵(CAR)-NKT cells or their derivatives collected from the stage 3 culture were resuspended in the expansion medium supplemented with 10 ng ml⁻¹ IL-7, 10 ng ml⁻¹ IL-15 and 1 μ g ml⁻¹ Ultra-LEAF Purified Anti-Human CD28 antibody (Clone CD28.2; BioLegend) at 5×10^5 cells per ml; 30 ml cell suspension was added into each plate. After 3 days of culture, cells were collected and resuspended in fresh expansion medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15, at $0.5\text{--}1 \times 10^6$ cells per ml. Cells were passaged 2–3 times per week to maintain a cell density at $0.5\text{--}1 \times 10^6$ cells per ml; fresh medium was added at every passage. As the cell population reaches a high number during the later stages, these cells can be transferred and cultured in G-Rex 6M Well Plates (Wilson Wolf).
- (2) The α GC-PBMC expansion approach. Healthy donor PBMCs were loaded with α -Galactosylceramide (α GC; Avanti Polar Lipids) at 5 μ g ml⁻¹ in C10 medium for 1 h following a previously established protocol²³. The resulting α GC-loaded PBMCs (α GC-PBMCs) were then irradiated at 6,000 rads using a Rad Source RS-2000 X-Ray Irradiator (Rad Source Technologies). Mature Allo¹⁵CAR-NKT cells and derivatives collected from the stage 3 culture were mixed with the irradiated α GC-PBMCs at a 1:5 ratio, resuspended in expansion medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15 at $0.5\text{--}1 \times 10^6$ cells per ml and seeded into the 150 mm cell culture dishes (Thermo Fisher Scientific) at 30 ml per plate. Cells were passaged 2–3 times per week to maintain a cell density at $0.5\text{--}1 \times 10^6$ cells per ml; fresh medium was added at every passage. As the cell population reaches a

high number during the later stages, these cells can be transferred and cultured in G-Rex 6M Well Plates (Wilson Wolf).

- (3) The aAPC expansion approach. aAPCs were irradiated at 10,000 rads using a Rad Source RS-2000 X-Ray Irradiator (Rad Source Technologies). Mature Allo¹⁵CAR-NKT cells and derivatives collected from the stage 3 culture were mixed with the irradiated aAPCs at a 1:1 ratio, resuspended in expansion medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15 at $0.5\text{--}1 \times 10^6$ cells per ml, and seeded into the 150 mm cell culture dishes (Thermo Fisher Scientific) at 30 ml per plate. Cells were passaged 2–3 times per week to maintain a cell density at $0.5\text{--}1 \times 10^6$ cells per ml; fresh medium was added at every passage. As the cell population reaches a high number, cells can be transferred to and cultured in G-Rex 6M Well Plates (Wilson Wolf).

Generation of PBMC-derived conventional α β T, NKT and NK cells
Healthy donor PBMCs were used to generate the PBMC-derived conventional α β T, NKT and NK cells (denoted as PBMC-Tc, PBMC-NKT and PBMC-NK cells, respectively).

To generate PBMC-Tc cells, PBMCs were stimulated with Dynabeads Human T-Activator CD3/CD28 (Thermo Fisher Scientific) according to the manufacturer's instructions, followed by culturing in the C10 medium supplemented with 20 ng ml⁻¹ IL-2 for 2–3 weeks.

To generate PBMC-NKT cells, PBMCs were sorted with magnetic-activated cell sorting (MACS) via Anti-iNKT Microbeads (Miltenyi Biotech) labeling to enrich NKT cells, following the manufacturer's instructions. The enriched NKT cells were mixed with donor-matched irradiated α GC-PBMCs at a ratio of 1:1, followed by culturing in C10 medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15 for 2–3 weeks. If needed, the resulting cultured cells could be further purified using fluorescence-activated cell sorting (FACS) via human NKT TCR antibody (Clone 6B11; BD Biosciences) staining.

To generate PBMC-NK cells, PBMCs were sorted with FACS using a FACS Aria III Sorter (BD Biosciences) via human CD56 antibody (Clone HCD56; BioLegend) labeling or with MACS using a Human NK Cell Isolation Kit (Miltenyi Biotech), following the manufacturers' instructions.

Generation of CAR-engineered conventional α β T (CAR-T) cells
Non-treated tissue culture 24-well plates (Corning) were coated with Ultra-LEAF Purified Anti-Human CD3 Antibody (Clone OKT3; BioLegend) at 1 μ g ml⁻¹ (500 μ l per well), at room temperature for 2 h or at 4 °C overnight. Healthy donor PBMCs were resuspended in the C10 medium supplemented with 1 μ g ml⁻¹ Ultra-LEAF Purified Anti-Human CD28 Antibody (Clone CD28.2, BioLegend) and 30 ng ml⁻¹ IL-2, followed by seeding in the precoated plates at 1×10^6 cells per ml (1 ml per well). On day 2, cells were transduced with lentiviruses for 24 h. The resulting CAR-T cells were expanded for about 2 weeks in C10 medium and cryopreserved for future use.

Generation of PBMC-derived CAR-engineered NKT (PBMC-CAR-NKT) cells

Healthy donor PBMCs were sorted with MACS via Anti-iNKT Microbeads (Miltenyi Biotech) labeling to enrich NKT cells, following the manufacturer's instructions. The enriched NKT cells were mixed with donor-matched irradiated α GC-PBMCs at a ratio of 1:1, followed by culturing in C10 medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15. On day 3, NKT cells were transduced with Lenti/BCAR-IL15 viruses for 24 h. The resulting CAR-NKT cells were expanded for about 2 weeks in C10 medium supplemented with 10 ng ml⁻¹ IL-7 and IL-15 and cryopreserved for future use.

In vitro tumor cell killing assay

Tumor cells (1×10^4 cells per well) were cocultured with therapeutic cells (at ratios indicated in figure legends) in Corning 96-well clear bottom black plates for 8–24 h, in C10 medium with or without the

addition of α GC (100 ng ml^{-1}). At the end of culture, live tumor cells were quantified by adding D-luciferin ($150 \mu\text{g ml}^{-1}$; Caliper Life Science) to cell cultures and reading out luciferase activities using an Infinite M1000 microplate reader (Tecan). In some experiments, $10 \mu\text{g ml}^{-1}$ LEAF purified antihuman NKG2D (Clone 1D11, BioLegend), antihuman DNAM-1 antibody (Clone 11A8, BioLegend) or LEAF purified mouse IgG2bk isotype control antibody (Clone MG2B-57, BioLegend) was added to cocultures, to study the NK activating receptor-mediated tumor cell killing mechanism.

In vitro assays using samples from patients with MM

Primary BM samples from patients with MM were collected and subsequently diluted in PBS and subjected to density gradient centrifugation using Ficoll-Paque (Thermo Fisher Scientific) to obtain mononuclear cells following the manufacturer's instructions. The resulting cells were cryopreserved for future use.

In one assay, the primary samples from patients with MM were analyzed for tumor cell phenotype and the TME composition using flow cytometry. Tumor cells were sorted using a Human Tumor Cell Isolation Kit (Miltenyi Biotec) and/or identified as $\text{CD45}^{\text{hi}}\text{CD31}^{\text{lo}}\text{FAP}$ (fibroblast activation protein)⁺ cells, T cells were identified as $\text{CD45}^{\text{hi}}\text{CD3}^{\text{hi}}$ cells, B cells were identified as $\text{CD45}^{\text{hi}}\text{CD19}^{\text{hi}}$ cells, NK cells were identified as $\text{CD45}^{\text{hi}}\text{CD56}^{\text{hi}}$ cells, monocytes and macrophages were identified as $\text{CD45}^{\text{hi}}\text{CD11b}^{\text{hi}}\text{CD14}^{\text{hi}}$ cells, granulocytes were identified as $\text{CD45}^{\text{hi}}\text{CD11b}^{\text{hi}}\text{CD15}^{\text{hi}}$ cells and granulocyte-like MDSCs were identified as $\text{CD45}^{\text{hi}}\text{CD15}^{\text{hi}}\text{HLA-DR}^{\text{lo}}$ cells. Surface expression of BCMA, CD1d and NK ligands on tumor or/and immune cells were also analyzed using flow cytometry.

In another assay, the primary samples from patients with MM were used to study tumor cell killing by $\text{Allo}^{15}\text{BCAR-NKT}$ cells. Tumor cells were presorted using a Human Tumor Cell Isolation Kit (Miltenyi Biotec), followed by coculturing with various therapeutic cells (therapeutic cell:tumor cell ratio 1:1) in C10 medium in Corning 96-well Round Bottom Cell Culture plates for 24 h. At the end of culture, cells were collected and live MM tumor cells (identified as $\text{CD45}^{\text{hi}}\text{CD3}^{\text{lo}}\text{6B11}^{\text{lo}}$) were analyzed using flow cytometry.

In another assay, the primary samples from patients with MM were used to study the TME targeting by $\text{Allo}^{15}\text{BCAR-NKT}$ cells. Patient samples were directly cocultured with $\text{Allo}^{15}\text{BCAR-NKT}$ cells (therapeutic cell:tumor cell ratio 1:1) in C10 medium in Corning 96-well Round Bottom Cell Culture plates for 24 h. At the end of culture, cells were collected and the TME targeting of $\text{Allo}^{15}\text{BCAR-NKT}$ cells was assessed using flow cytometry by quantifying live human monocytes and macrophages (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{hi}}\text{CD14}^{\text{hi}}\text{CD11b}^{\text{lo}}$), MDSCs (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{med}}\text{CD15}^{\text{hi}}\text{CD11b}^{\text{lo}}\text{HLA-DR}^{\text{lo}}$), CD4 T cells (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{hi}}\text{CD3}^{\text{hi}}\text{CD4}^{\text{lo}}$), CD8 T cells (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{hi}}\text{CD3}^{\text{hi}}\text{CD8}^{\text{hi}}$), B cells (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{hi}}\text{CD3}^{\text{lo}}\text{CD19}^{\text{hi}}$), granulocytes (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{med}}\text{CD15}^{\text{hi}}\text{CD11b}^{\text{lo}}\text{HLA-DR}^{\text{lo}}$) and NK cells (identified as $6\text{B11}^{\text{lo}}\text{CD45}^{\text{hi}}\text{CD3}^{\text{lo}}\text{CD56}^{\text{hi}}$).

In vitro MLR assay for studying GvH response

PBMCs of multiple random healthy donors were irradiated at 2,500 rads and used as stimulators to study the GvH response of $\text{Allo}^{15}\text{BCAR-NKT}$ cells as responders. BCAR-T cells were included as a responder control. Stimulators (5×10^5 cells per well) and responders (2×10^4 cells per well) were cocultured in 96-well round-bottom plates in C10 medium for 4 days; the cell culture supernatants were then collected to measure IFN γ production using enzyme-linked immunosorbent assay (ELISA).

In vitro MLR assay for studying T cell-mediated allorejection

PBMCs of multiple healthy donors were used as responders to study the T cell-mediated allorejection of $\text{Allo}^{15}\text{BCAR-NKT}$ cells as stimulators (irradiated at 2,500 rads). PBMC-derived BCAR-T cells were included as a stimulator control. Stimulators (5×10^5 cells per well) and responders

(2×10^4 cells per well) were cocultured in 96-well round-bottom plates in C10 medium for 4 days; the cell culture supernatants were then collected to measure IFN γ production using ELISA.

In vitro MLR assay for studying NK cell-mediated allorejection

PBMC-derived NK cells obtained from multiple healthy donors were employed to investigate the NK cell-mediated allorejection of $\text{Allo}^{15}\text{BCAR-NKT}$ cells. Conventional BCAR-T cells were included as an allogeneic subject control. PBMC-NK cells (2×10^4 cells per well) and the corresponding allogeneic subject cells (2×10^4 cells per well) were cocultured in 96-well round-bottom plates with C10 medium for 24 h. Subsequently, the cell cultures were collected to quantify live cells via flow cytometry.

In vivo BLI

Bioluminescence live animal imaging (BLI) was performed using a Spectral Advanced Molecular Imaging (AMI) HTX imaging system (Spectral Instrument Imaging). Live animal images were acquired 5 min after intraperitoneal (i.p.) injection of D-luciferin for total body bioluminescence. To monitor the signal of tumor cells, 1 mg per mouse of D-luciferin was injected. To monitor the signal of therapeutic cells, 3 mg per mouse of D-luciferin was injected. Imaging data were analyzed using AURA imaging software (Spectral Instrument Imaging, v.3.2.0).

In vivo antitumor efficacy study of $\text{Allo}^{15}\text{BCAR-NKT}$ cells: human MM xenograft NSG mouse model

Experimental design is shown in Fig. 3a. Briefly, on Day 0, NSG mice received intravenously (i.v.) inoculation of MM-FG human MM cells (1×10^6 cells per mouse). On day 10, the experimental mice received i.v. injection of vehicle (100 μl PBS per mouse), $\text{Allo}^{15}\text{BCAR-NKT}$ cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse), or control BCAR-T cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse). Over the experiment, mice were monitored for survival and their tumor loads were measured using BLI.

In vivo antitumor efficacy study of $\text{Allo}^{15}\text{CAR19-NKT}$ cells: human B cell leukemia xenograft NSG mouse model

Experimental design is shown in Supplementary Fig. 15h. Briefly, on day 0, NSG mice received i.v. inoculation of Nalm6-FG human B cell leukemia cells (1×10^6 cells per mouse). On day 4, the experimental mice received an i.v. injection of vehicle (100 μl PBS per mouse), $\text{Allo}^{15}\text{CAR19-NKT}$ cells (5×10^6 CAR⁺ cells in 100 μl PBS per mouse) or control CAR19-T cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse). Over the experiment, mice were monitored for survival and their tumor loads were measured using BLI.

In vivo PK-PD study of $\text{Allo}^{15}\text{BCAR-NKT/FG}$ cells: human MM xenograft NSG mouse model

Experimental design is shown in Fig. 3f. Briefly, on day 0, NSG mice received i.v. inoculation of MM.1S cells (1×10^6 cells per mouse; denoted as MM). On day 10, the experimental mice received i.v. injection of $\text{Allo}^{15}\text{BCAR-NKT/FG}$ cells (10×10^6 CAR⁺ cells in 100 μl of PBS per mouse) or control BCAR-T/FG cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse). Over the experiment, mice were monitored for survival and their therapeutic cells were measured using BLI. Note in this study, therapeutic immune cells, but not the tumor cells, were labeled with FG.

In vivo CRS study

Experimental design is shown in Fig. 5p. Briefly, on day 0, NSG mice received i.v. inoculation of MM-FG cells (5×10^6 cells per mouse). On day 10, the experimental mice received i.v. injection of vehicle (100 μl PBS per mouse), $\text{Allo}^{15}\text{BCAR-NKT}$ cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse) or control BCAR-T cells (10×10^6 CAR⁺ cells in 100 μl PBS per mouse). On days 11 and 13, blood samples were collected from the experimental mice, and their serum IL-6 and SAA-3 were measured

using ELISA. A Mouse SAA-3 ELISA Kit (Millipore Sigma) was used to measure SAA-3, following the manufacturer's instructions.

Single-cell TCR sequencing

PBMC-Tc, PBMC-NKT, $Allo^0$ NKT and $Allo^{15}$ BCAR-NKT cells were sorted using a FACSAria III Sorter (BD Biosciences). Sorted cells were immediately delivered to the UCLA Technology Center for Genomics and Bioinformatics (TCGB) Core to perform single-cell TCR sequencing using a 10X Genomics Chromium Controller Single Cell Sequencing System (10X Genomics), following the manufacturer's instructions and the TCGB Core's standard protocols. Libraries were constructed using an Illumina TruSeq RNA Sample Prep Kit (Illumina) and sequenced with 150 bp paired-end reads (5,000 reads per cell) on an Illumina NovaSeq 6000 sequencer. The reads were mapped to the human TCR reference genome (hg38) using Cell Ranger VDJ (10X Genomics). The frequencies of the α or β chain recombination were plotted.

Bulk RNA-seq

A total of 24 cell samples were analyzed, including three $Allo^0$ NKT, three $Allo^0$ BCAR-NKT, three $Allo^{15}$ BCAR-NKT, two $Allo^0$ CAR19-NKT, three PBMC-NKT, eight PBMC-Tc and two PBMC-NK cell samples; the numbers indicate different donors. Note that $Allo^{15}$ (CAR)-NKT cells displayed a $CD4^+CD8^{+/-}$ phenotype, therefore the $CD4^+$ populations of PBMC-NKT and PBMC-Tc cells were used in this study.

Cell samples were sorted using a FACSAria III Sorter (BD Biosciences). Total RNAs were isolated from each cell sample using a miRNeasy Mini Kit (Qiagen) and delivered to the UCLA TCGB Core to perform bulk RNA-seq using an Illumina HiSeq3000, following the manufacturer's instructions and the TCGB Core's standard protocols. Complementary DNAs were synthesized using an iScript cDNA Synthesis Kit (BioRad). Libraries were constructed using an Illumina TruSeq Stranded Total RNA Sample Prep kit and sequenced with 50 bp single-end reads (targeting 20×10^6 reads per sample) on an Illumina HiSeq3000. The raw sequencing data underwent quality control and adapter trimming using fastp (v.0.23.2), and trimmed reads were mapped to the human genome (hg38) using with STAR v.2.7.9a with default parameters. The read counts for each gene were obtained using featureCounts from Subread package (v.2.0.3). Sequencing depth normalized counts (in counts per million) were obtained using edgeR (v.3.38.4), and the log transformation with psuedocount 1 was used for principal component analysis.

scRNA-seq

In one study, scRNA-seq was used to investigate the lineage commitment of $Allo^{15}$ BCAR-NKT cells during their in vitro differentiation from gene-engineered HSPCs. Cord-blood $CD34^+$ HSPCs of the same donor were transduced with either a Lenti/iNKT-BCAR or Lenti/iNKT-BCAR-IL15 lentivector, followed by differentiating into $Allo^0$ BCAR-NKT or $Allo^{15}$ BCAR-NKT cell product, respectively, in the Ex Vivo HSPC-Derived CAR-NKT Cell Culture. At days 7, 14, 21, 28 and 42, cell cultures were collected for scRNA-seq.

In another study, scRNA-seq was used to examine the in vivo gene profiles of therapeutic cells, tumor cells and mouse immune cells in an MM-FG human MM xenograft NSG mouse model as shown in Fig. 4a. The three types of therapeutic cell ($Allo^{15}$ BCAR-NKT, $Allo^0$ BCAR-NKT, and BCAR-T cells) were collected at three time points (days 0, 28 and 35). Day 0 samples were collected before infusion, day 28 samples were collected from the BM tumor sites of the experimental mice and day 35 samples were generated by stimulating day 28 samples in vitro with BCMA-expressing aAPCs for 1 week. The tumor cells (MM-FG) and mouse immune cells (mouse $CD45^+$ cells) were collected on day 28. Cells isolated from ten mice of each experimental group were combined for scRNA-seq.

Freshly collected samples were immediately delivered to the UCLA TCGB Core for library construction and scRNA-seq. Cells were

quantified using a Cell Countess II automated cell counter (Invitrogen/Thermo Fisher Scientific). A total of 10,000 cells from each experimental group were loaded on the Chromium platform (10X Genomics), and libraries were constructed using the Chromium Next GEM Single Cell 3' Kit v3.1 and the Chromium Next GEM Chip G Single Cell Kit (10X Genomics), according to the manufacturer's instructions. Library quality was assessed using the D1000 ScreenTape on a 4200 TapeStation System (Agilent Technologies). Libraries were sequenced on an Illumina NovaSeq using the NovaSeq S4 Reagent Kit (100 cycles; Illumina).

For cell clustering and annotation, the merged digital expression matrix generated by Cell Ranger was analyzed using an R package Seurat (v.4.0.0) following the guidelines. Briefly, after filtering the low-quality cells, the expression matrix was normalized using NormalizeData function, followed by selecting variable features across datasets using FindVariableFeatures and SelectIntegrationFeatures functions. To correct the batch effect, FindIntegrationAnchors and IntegrateData functions were used based on the selected feature genes. The corrected dataset was subjected to standard Seurat workflow for dimension reduction and clustering. In this study, clusters of therapeutic cells were manually merged and annotated based on gene signatures reported from Human Protein Atlas (proteinatlas.org) and previous studies³⁸⁻⁴⁵, and clusters of mouse immune cells were merged and annotated based on the immune lineage markers. AddModuleScore was used to calculate module scores of each list of gene signatures, and FeaturePlot function was used to visualize the expression of each signature in the UMAP plots.

For gene set enrichment analysis (GSEA), clusterProfiler packages^{87,88} were used to calculate the enrichment scores of each cluster in the signature gene list.

For RNA-velocity analysis, the .LOOM files containing spliced and unspliced expression matrices were generated for each sample. The further analysis was conducted using the Velocity.R-package (v.0.6)⁸⁹. After loading .LOOM file information through the ReadVelocity function, databases were merged and the RunVelocity function was executed to obtain the velocity vectors. Finally, the velocities were projected into a lower-dimensional embedding using the velocity_graph function and visualized on the UMAP embedding in each intended cell cluster using the show.velocity.on.embedding.cor function. All velocity functions were used with default parameters.

Methylation sequencing

Genomic DNA was isolated from experimental samples using a QIA-GEN DNeasy Blood & Tissue kit, then sonicated using a Covaris M220 Focused-ultrasonicator, following the manufacturers' instructions. DNA fragments of around 250 bp were enriched using the Ampure XP beads (Beckman-Coulter), then subjected to DNA library preparation using an NEBNext Ultra II DNA library prep Kit (cat. no. E7645) following the manufacturers' instructions. The DNA libraries were then subjected to sequencing on Illumina NovaSeq sequencer with 2×150 bp configuration (Azenta). The raw bisulfite sequencing reads were trimmed by cutadapt to remove sequencing adapters. The trimmed reads were then aligned to hg19 reference genome by Bismark⁹⁰. Then duplicated reads from PCR amplification were identified and removed by Bismark. The deduplicated reads were then sorted and indexed using Samtools. After that, the methylated and unmethylated cytosines were counted at every CpG site by Bismark. The methylation at a gene promoter region was quantified as the beta value, that is, the ratio between the number of methylated cytosines and the total number of cytosines mapped to the region. The promoter region of a gene was identified by GeneHancer with the highest confidence score. The beta values were then used for heatmap visualization.

Histology

Tissues (that is, heart, liver and lung) were collected from experimental mice, fixed in 10% Neutral Buffered Formalin for up to 36 h and

embedded in paraffin for sectioning (5 μm thickness). Tissue sections were prepared and stained with hematoxylin and eosin by the UCLA Translational Pathology Core Laboratory, following the Core's standard protocols. Stained sections were imaged using an Olympus BX51 upright microscope equipped with an Optronics Macrofire CCD camera (AU Optronics). The images were analyzed using Optronics PictureFrame software (AU Optronics).

Western blot

Western blot was used to analyze the IFN γ signaling events in Allo^{15} BCAR-NKT cells. Conventional BCAR-T cells and quiescent T cells prepared from the same healthy donor PBMCs were included as controls. Quiescent T cells were sorted with MACS from PBMCs using mixed antihuman CD4 and CD8 magnetic beads (Miltenyi Biotec) to avoid the potential CD3 stimulation, following the manufacturer's instructions.

Cells were stimulated with IFN γ (10 ng ml $^{-1}$) for 15 min (for blotting p-STAT1/STAT1), 18 h (for blotting IRF-1) or 48 h (for blotting CIITA/SPI). Total proteins were extracted using a RIPA lysis buffer (Thermo Fisher Scientific) containing 20 mM HEPES (pH 7.6), 150 mM NaCl, 1 mM EDTA, 1% TritonX-100 and protease-phosphatase inhibitor cocktail (Cell Signaling Technology). Protein concentration was measured using a Bicinchoninic Acid Assay Kit (Thermo Fisher Scientific). Equal amounts of total protein were resolved on a 4–15% Mini-PROTEAN TGX Precast Protein Gel (BioRad) and then transferred to a polyvinylidene difluoride membrane by electrophoresis. The following antibodies were used to blot for the proteins of interest: antihuman CIITA (Invitrogen, cat. no. PA5-21031, 1:1,000 dilution), antihuman p-STAT1(Y701) (Clone 58D6, Cell Signaling Technology, CST, cat. no. 9167, 1:1,000 dilution), antihuman STAT1 (Clone DIK9Y, CST, cat. no. 14994, 1:1,000 dilution), antihuman IRF-1 (Clone D5E4, CST, cat. no. 8478, 1:1,000 dilution), antihuman SPI (Clone D4C3, CST, cat. no. 9389, 1:1,000 dilution) and secondary antirabbit IgG (CST, cat. no. 7074, 1:1,000 dilution). β -Actin (Clone D6A8, CST, cat. no. 8457, 1:3,000 dilution) and glyceraldehyde 3-phosphate dehydrogenase (Clone D16H11, CST, cat. no. 5174, 1:3,000 dilution) were used as internal controls. Signals were visualized using a ChemiDoc Imaging Systems (Bio-Rad). The data were analyzed using ImageJ (v.1.53s).

Generation and analysis of Allo^{15} BCAR-NKT cells overexpressing a constitutively active STAT1 (denoted as Allo^{15} BCAR-NKT/STATIC cells)

Allo^{15} BCAR-NKT/STATIC cells were generated by further engineering the Allo^{15} BCAR-NKT cells to overexpress a constitutively active human STAT1 (STATIC)⁸⁵. The Lenti/STATIC-EGFP vector or a mock control Lenti/Fluc-EGFP vector was added into the Allo^{15} BCAR-NKT cell culture during stage 4 (on day 3 of stage 4), followed by the continuation and completion of stage 4 culture. The resulting Allo^{15} BCAR-NKT/STATIC cells, as well as their mock control Allo^{15} BCAR-NKT/FG cells, were then analyzed using flow cytometry. The overexpression of transgenic STATIC was measured by EGFP coexpression in Allo^{15} BCAR-NKT/STATIC cells. The rescue effect of STATIC overexpression in EGFP $^{+}$ Allo^{15} BCAR-NKT/STATIC cells was assessed by comparing these cells with the control EGFP $^{+}$ Allo^{15} BCAR-NKT/FG cells for cell surface expression of HLA-I/II molecules and NK ligands.

Statistics

Graphpad Prism v.8 software (Graphpad) was used for statistical data analysis. Student's two-tailed *t*-test was used for pairwise comparisons. Ordinary one-way analysis of variance (ANOVA) followed by Tukey's or Dunnett's multiple comparisons test was used for multiple comparisons. A log rank (Mantel-Cox) test adjusted for multiple comparisons was used for Meier survival curves analysis. Data are presented as the mean \pm s.e.m., unless otherwise indicated. In all figures and figure legends, *n* represents the number of samples or animals used in the indicated experiments. A *P* value of less than 0.05 was considered significant and NS denotes not significant.

Reporting summary

Further information on research design is available in the Nature Portfolio Reporting Summary linked to this article.

Data availability

All data associated with this study are present in the paper or Supplementary Information. The genomics data generated during this study are available from the public repository Gene Expression Omnibus Database: GSE253982 (scTCR-seq, related to Fig. 1h), GSE245375 (bulk RNA-seq, related to Fig. 1i), GSE241996 (scRNA-seq for the in vitro study related to Fig. 1g), GSE241998 (scRNA-seq for the in vivo study related to Fig. 4), GSE241999 (scRNA-seq for the in vivo study related to Fig. 5a–c), GSE241997 (scRNA-seq for the in vivo study related to Fig. 5h–n) and GSE239648 (Methyl-seq, related to Fig. 6n). Additional information and materials will be made available upon reasonable request. Source data are provided with this paper.

References

- Wu, C. et al. STAT1 is phosphorylated and downregulated by the oncogenic tyrosine kinase NPM-ALK in ALK-positive anaplastic large-cell lymphoma. *Blood* **126**, 336–345 (2015).
- Seet, C. S. et al. Generation of mature T cells from human hematopoietic stem and progenitor cells in artificial thymic organoids. *Nat. Methods* **14**, 521–530 (2017).
- Yu, G., Wang, L.-G., Han, Y. & He, Q.-Y. clusterProfiler: an R package for comparing biological themes among gene clusters. *OMICS* **16**, 284–287 (2012).
- Wu, T. et al. clusterProfiler 4.0: a universal enrichment tool for interpreting omics data. *Innovation* **2**, 100141 (2021).
- La Manno, G. et al. RNA velocity of single cells. *Nature* **560**, 494–498 (2018).
- Krueger, F. & Andrews, S. R. Bismark: a flexible aligner and methylation caller for Bisulfite-Seq applications. *Bioinformatics* **27**, 1571–1572 (2011).

Acknowledgements

We thank the UCLA animal facility for providing animal support; the UCLA Translational Pathology Core Laboratory for providing histology support; the UCLA TCGB facility for providing RNA-seq services; the UCLA CFAR Virology Core for providing human cells; and the UCLA Broad Stem Cell Research Centre (BSCRC) Flow Cytometry Core Facility for cell sorting support. This work was supported by a Partnering Opportunity for Discovery Stage Research Projects Award and a Partnering Opportunity for Translational Research Projects Awards from the California Institute for Regenerative Medicine (grant nos. DISC2-11157 and TRAN1-12250, to L.Y.), a Department of Defense CDMRP PRCRP Impact Award (grant no. CA200456 to L.Y.), a UCLA BSCRC Innovation Award (to L.Y.), and an Ablon Scholars Award (to L.Y.). Y.-R.L. is a postdoctoral fellow supported by a UCLA MIMG M. John Pickett Post-Doctoral Fellow Award and a CIRM-BSCRC Postdoctoral Fellowship. Yang Zhou is a predoctoral fellow supported by the UCLA Dissertation Year Fellowship. J.Y. is a predoctoral fellow supported by the UCLA BSCRC Predoctoral Fellowship. D.L. and Z.L. are postdoctoral fellows supported by T32 Tumor Immunology Training Grant postdoctoral fellowships (USHHS Ruth L. Kirschstein Institutional National Research Service Award, grant no. T32-CA009120).

Author contributions

Y.-R.L., Yang Zhou, J.Y. and L.Y. designed the experiments, analyzed the data and wrote the paper. L.Y. conceived and oversaw the study, with assistance from Y.-R.L., Yang Zhou and J.Y., and suggestions from D.B.K., P.W., X.J.Z. and M.P. Y.-R.L., Yang Zhou and J.Y. performed all experiments, with assistance from Y.J.K., M.L., D.L., K.Z., Y.C., Y. Zhu, Y.-C.W., Z.L., Y.Y., Z.S.D., X.C., T.H., A.B., A.K., M.W., Y.F., J.H., Y. Zhang

and Z.H. Y.-C.W., W.G., F.M., C.P. and A.J.L. helped with analysis of RNA-seq data. S.L. and Yonggang Zhou helped with the study of DNA methylation. E.Z. helped with SEM imaging. D.B.K. provided guidance on development of lentiviral vectors. J.J.Z. helped with the statistical analysis of data. C.S.S. provided the aAPC line. B.R.P. and S.M.L. provided and assisted on the analysis of primary BM samples from patients with MM.

Competing interests

Y.-R.L., Yang Zhou, J.Y., Y.J.K., D.L., Z.L., P.W. and L.Y. are inventors on patents relating to this study filed by UCLA. Yang Zhou is currently an employee of Amberstone Biosciences. J.Y. is currently an employee of Appia Bio. Y.J.K. is currently an employee of Nkarta. Z.L. is currently an employee of Allogene. Y.Y. is currently an employee of Retro Biosciences. Z.S.D. is currently an employee of Guggenheim Securities. F.M. is currently an employee of Amgen. C.S.S. is a cofounder and stockholder of Pluto Immunotherapeutics. D.B.K. is a stockholder and scientific advisory board member of Allogene Therapeutics, Myogene Bio, ImmunoVec and Pluto Immunotherapeutics. P.W. is a cofounder, stockholder and advisory board member of Simnova Bio, TCRCure Biopharma, Appia Bio, and is

a scientific advisor to Grit Biotechnolgy. L.Y. is a scientific advisor to AlzChem and Amberstone Biosciences, and a cofounder, stockholder and advisory board member of Appia Bio. Appia Bio licensed some patents relating to this study from UCLA. None of the declared companies contributed to or directed any of the research reported in this article. The other authors declare no competing interests.

Additional information

Supplementary information The online version contains supplementary material available at <https://doi.org/10.1038/s41587-024-02226-y>.

Correspondence and requests for materials should be addressed to Lili Yang.

Peer review information *Nature Biotechnology* thanks Paolo Dellabona, Leonid Metelitsa and the other, anonymous, reviewer(s) for their contribution to the peer review of this work.

Reprints and permissions information is available at www.nature.com/reprints.

Reporting Summary

Nature Portfolio wishes to improve the reproducibility of the work that we publish. This form provides structure for consistency and transparency in reporting. For further information on Nature Portfolio policies, see our [Editorial Policies](#) and the [Editorial Policy Checklist](#).

Statistics

For all statistical analyses, confirm that the following items are present in the figure legend, table legend, main text, or Methods section.

- | | |
|-------------------------------------|--|
| n/a | Confirmed |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> The exact sample size (n) for each experimental group/condition, given as a discrete number and unit of measurement |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A statement on whether measurements were taken from distinct samples or whether the same sample was measured repeatedly |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> The statistical test(s) used AND whether they are one- or two-sided
<i>Only common tests should be described solely by name; describe more complex techniques in the Methods section.</i> |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A description of all covariates tested |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A description of any assumptions or corrections, such as tests of normality and adjustment for multiple comparisons |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> A full description of the statistical parameters including central tendency (e.g. means) or other basic estimates (e.g. regression coefficient) AND variation (e.g. standard deviation) or associated estimates of uncertainty (e.g. confidence intervals) |
| <input type="checkbox"/> | <input checked="" type="checkbox"/> For null hypothesis testing, the test statistic (e.g. F , t , r) with confidence intervals, effect sizes, degrees of freedom and P value noted
<i>Give P values as exact values whenever suitable.</i> |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> For Bayesian analysis, information on the choice of priors and Markov chain Monte Carlo settings |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> For hierarchical and complex designs, identification of the appropriate level for tests and full reporting of outcomes |
| <input checked="" type="checkbox"/> | <input type="checkbox"/> Estimates of effect sizes (e.g. Cohen's d , Pearson's r), indicating how they were calculated |

Our web collection on [statistics for biologists](#) contains articles on many of the points above.

Software and code

Policy information about [availability of computer code](#)

Data collection: MACSQuant Analyzer 10 Flow Cytometer, Illumina NovaSeq 6000 sequencer, Infinite M1000 microplate reader (Tecan), Spectral Advanced Molecular Imaging (AMI) HTX imaging system (Spectral instrument Imaging), ChemiDoc™ Imaging Systems (BIO-RAD)

Data analysis: FlowJo software (version 9), AURA imaging software (version 3.2.0), Graphpad Prism (version 8), Image J (version 1.53s), R package Seurat (v.4.0.0)

For manuscripts utilizing custom algorithms or software that are central to the research but not yet described in published literature, software must be made available to editors and reviewers. We strongly encourage code deposition in a community repository (e.g. GitHub). See the Nature Portfolio [guidelines for submitting code & software](#) for further information.

Data

Policy information about [availability of data](#)

All manuscripts must include a [data availability statement](#). This statement should provide the following information, where applicable:

- Accession codes, unique identifiers, or web links for publicly available datasets
- A description of any restrictions on data availability
- For clinical datasets or third party data, please ensure that the statement adheres to our [policy](#)

All data associated with this study are present in the paper or Supplemental Information. The genomics data generated during this study are available from the public repository Gene Expression Omnibus Database: GSE (scTCR-seq, related to Fig. 1h), GSE245375 (Bulk RNA-seq, related to Fig. 1i), GSE241996 (scRNA-seq for

the in vitro study related to Fig. 1g), GSE241998 (scRNA-seq for the in vivo study related to Fig. 4), GSE241999 (scRNA-seq for the in vivo study related to Fig. 5a-5c), GSE241997 (scRNA-seq for the in vivo study related to Fig. 5h-5n), GSE239648 (Methy-seq, related to Fig. 6n). Source data are provided with this paper. Additional information and materials will be made available upon reasonable request.

Human research participants

Policy information about [studies involving human research participants and Sex and Gender in Research](#).

Reporting on sex and gender	Sex and gender are not considered in the design of this study.
Population characteristics	<p>Patient #1 Collection data: 2/13/2023 Age: 64 Diagnosis: Multiple myeloma FISH: gain of 1q, gain of 5p, monosomy 13, and loss of 17p (high-risk cytogenetics) Number of prior lines of therapy: 6</p> <p>Patient #2 Collection data: 3/06/2023 Age: 65 Diagnosis: Multiple myeloma FISH: gain of 1q and CCND1-IGH fusions [t(11;14)]. Number of prior lines of therapy: 4</p> <p>Patient #3 Collection data: 3/29/2023 Age: 61 Diagnosis: Multiple myeloma FISH: +5, +11q, and 16q-. Number of prior lines of therapy: 5</p>
Recruitment	Healthy donor PBMCs were obtained from the UCLA/CFAR Virology Core Laboratory. Primary MM patient bone marrow samples were collected at the Ronald Reagan UCLA Medical Center from consented patients and processed.
Ethics oversight	Healthy donor PBMCs were obtained from the UCLA/CFAR Virology Core Laboratory without identification information under federal and state regulations. Primary MM patient bone marrow samples were collected at the Ronald Reagan UCLA Medical Center from consented patients through an IRB-approved protocol (IRB#21-001444).

Note that full information on the approval of the study protocol must also be provided in the manuscript.

Field-specific reporting

Please select the one below that is the best fit for your research. If you are not sure, read the appropriate sections before making your selection.

Life sciences Behavioural & social sciences Ecological, evolutionary & environmental sciences

For a reference copy of the document with all sections, see [nature.com/documents/nr-reporting-summary-flat.pdf](https://www.nature.com/documents/nr-reporting-summary-flat.pdf)

Life sciences study design

All studies must disclose on these points even when the disclosure is negative.

Sample size	For all studies samples size were defined on the basis of previous studies in the laboratory (PMID: 34841295, 37938576, and 34112755). For in vivo studies, between 3 to 10 mice per group were used dependent on the experimental design and were able to reach a significant level of 0.05. For in vitro studies, at least 3 biological group were used as cell donors with technical replicates. For studies involving primary MM patient samples, at least 3 patient samples per group were used.
Data exclusions	No data were excluded from the analyses.
Replication	All experiments were performed independently. Replicates of each individual experiment was stated in its figure legends. The experimental findings were reproduced with similar results. For single cell RNA sequencing, a mixture of 10 experimental mice were combined for each group.
Randomization	For in vivo studies, all mice received tumor cells at Day 0 then were randomized to split into different groups prior to therapeutic cells injection. All the in vitro and in vivo experiments each group was selected randomly.
Blinding	Experiments were not performed in a blinded fashion. All data were analyzed by software with objective standard.

Reporting for specific materials, systems and methods

We require information from authors about some types of materials, experimental systems and methods used in many studies. Here, indicate whether each material, system or method listed is relevant to your study. If you are not sure if a list item applies to your research, read the appropriate section before selecting a response.

Materials & experimental systems

n/a	Involved in the study
<input type="checkbox"/>	<input checked="" type="checkbox"/> Antibodies
<input type="checkbox"/>	<input checked="" type="checkbox"/> Eukaryotic cell lines
<input checked="" type="checkbox"/>	<input type="checkbox"/> Palaeontology and archaeology
<input type="checkbox"/>	<input checked="" type="checkbox"/> Animals and other organisms
<input checked="" type="checkbox"/>	<input type="checkbox"/> Clinical data
<input checked="" type="checkbox"/>	<input type="checkbox"/> Dual use research of concern

Methods

n/a	Involved in the study
<input checked="" type="checkbox"/>	<input type="checkbox"/> ChIP-seq
<input type="checkbox"/>	<input checked="" type="checkbox"/> Flow cytometry
<input checked="" type="checkbox"/>	<input type="checkbox"/> MRI-based neuroimaging

Antibodies

Antibodies used

Fluorochrome-conjugated antibodies specific for human CD45 (Clone H130, CAT#304026, 1:500 dilution), TCR $\alpha\beta$ (Clone I26, CAT#306716, 1:25 dilution), CD3 (Clone HIT3a, CAT#300329, 1:500 dilution), CD4 (Clone OKT4, CAT#317414, 1:400 dilution), CD8 (Clone SK1, CAT#344714, 1:500 dilution), CD45RO (Clone UCHL1, CAT#304216, 1:200 dilution), CD45RA (Clone HI100, CAT#304105, 1:5000 dilution), CD161 (Clone HP-3G10, CAT#339928, 1:50 dilution), CD69 (Clone FN50, CAT#310909, 1:50 dilution), CD56 (Clone HCD56, CAT#362545, 1:10 dilution), CD1d (Clone 51.1, CAT#350308, 1:50 dilution), BCMA (19F2, CAT#357503, 1:50 dilution), CD14 (Clone HCD14, CAT#325608, 1:100 dilution), CD19 (Clone HIB19, CAT#363005, 1:100 dilution), CD11b (Clone ICRF44, CAT#301330, 1:500 dilution), CD15 (W6D3, CAT#323021, 1:500 dilution), CD112 (Clone TX31, CAT#337409, 1:200 dilution), CD155 (Clone SKII.4, CAT#337613, 1:200 dilution), MICA/MICB (Clone 6D4, CAT#320908, 1:50 dilution), Ganglioside GD2 (14G2a, CAT#357323, 1:50 dilution), NKG2D (Clone 1D11, CAT#320812, 1:50 dilution), DNAM-1 (Clone 11A8, CAT#338312, 1:50 dilution), CD158 (KIR2DL1/S1/S3/S5) (Clone HP-MA4, CAT#339510, 1:50 dilution), NKp30 (Clone P30-15, CAT#325207, 1:50 dilution), NKp44 (Clone P44-8, CAT#325107, 1:50 dilution), CD16 (Clone 3G8, CAT#302011, 1:50 dilution), NKG2A (Clone S19004C, CAT#375103, 1:50 dilution), NKG2C (Clone S19005E, CAT#375003, 1:50 dilution), CD62L (Clone P44-8, CAT#304813, 1:50 dilution), CD134 (Clone Ber-ACT35, CAT#350008, 1:50 dilution), LEF1 (Clone W17021C, CAT#621051, 1:50 dilution), CD8A (Clone C8/144B, CAT#372092, 1:500 dilution), CD8B (Clone QA20A40, CAT#387305, 1:500 dilution), IFN- γ (Clone B27, CAT#506518, 1:50 dilution), Granzyme B (Clone QA16A02, CAT#372204, 1:4000 dilution), Perforin (Clone dG9, CAT#308126, 1:50 dilution), TNF- α (Clone Mab11, CAT#502912, 1:4000 dilution), IL-2 (Clone MQ1-17H12, CAT#500341, 1:200 dilution), β 2-microglobulin (B2M) (Clone 2M2, CAT#316312, 1:5000 dilution), HLA-DR (Clone L243, CAT#307618, 1:250 dilution), and HLA-DR, DP, DQ (Clone Tü 39, CAT#361707, 1:250 dilution) were purchased from BioLegend. Fluorochrome-conjugated antibodies specific for mouse CD45 (Clone S18009F, CAT#157607, 1:5000 dilution), GR1 (Clone RB6-8C5, CAT#108411, 1:1000 dilution), CD1d (Clone 1B1, CAT#123521, 1:100 dilution), CD11b (Clone M1/70, CAT#101205, 1:5000 dilution), Sca1 (Clone D7, CAT#108111, 1:100 dilution), FIt3 (Clone A2F10, CAT#135305, 1:100 dilution), SLAM (Clone TC15-12F12.2, CAT#115913, 1:50 dilution), c-Kit (Clone 2B8, CAT#105825, 1:50 dilution) were purchased from BioLegend. Fluorochrome-conjugated antibody specific for human Glypican 3 (GPC3; Clone 024, CAT#ab275695, 1:200 dilution) was purchased from Abcam. Fluorochrome-conjugated antibodies specific for human CD34 (Clone 581, CAT#555822, 1:100 dilution) and human iNKT TCR Va24-J β 18 (Clone 6B11, CAT#552825, 1:10 dilution) were purchased from BD Biosciences. Fluorochrome-conjugated antibody specific for human iNKT TCR V β 11 (Clone C21, CAT#A66905, 1:50 dilution) was purchased from Beckman-Coulter. Fluorochrome-conjugated antibodies specific for human ULBP-1 (Clone 170818, CAT#FAB1380P, 1:50 dilution) and ULBP-2,5,6 (Clone 165903, CAT#FAB1298A5, 1:50 dilution) were purchased from R&D Systems. A goat anti-mouse IgG F(ab')₂ secondary antibody (CAT#31803, 1:50 dilution) was purchased from ThermoFisher. Fixable Viability Dye eFluor506 (e506, CAT#65-0866-14, 1:500 dilution) was purchased from Affymetrix eBioscience; mouse Fc Block (anti-mouse CD16/32, CAT#553142, 1:50 dilution) was purchased from BD Biosciences; and human Fc Receptor Blocking Solution (TrueStain FcX, CAT#422302, 1:100 dilution) was purchased from BioLegend.

The following antibodies were used to blot for the proteins of interest: anti-human CIITA (Invitrogen, CAT#PA5-21031, 1:1000 dilution), anti-human p-STAT1(Y701) (Clone 58D6, Cell Signaling Technology, CST, CAT#9167, 1:1000 dilution), anti-human STAT1 (Clone D1K9Y, CST, CAT#14994, 1:1000 dilution), anti-human IRF-1 (Clone D5E4, CST, CAT#8478, 1:1000 dilution), anti-human SP1 (Clone D4C3, CST, CAT#9389, 1:1000 dilution), and secondary anti-rabbit IgG (CST, CAT#7074, 1:1000 dilution). Beta-Actin (Clone D6A8, CST, CAT#8457, 1:3000 dilution) and GAPDH (Clone D16H11, CST, CAT#5174, 1:3000 dilution) were used as internal controls.

Mouse IL-6 was quantified by ELISA with paired purified anti-mouse IL-6 antibody (CAT#504501, 1:250 dilution) and biotin anti-mouse IL-6 antibody (CAT#504601, 1:500 dilution). These antibodies were purchased from BioLegend.

Validation

All the antibodies used are from commercial sources and have been validated by the vendors. Validation data are available on the manufacturer's website. Each antibody purchased for this study has relevant citation and profile.

Eukaryotic cell lines

Policy information about [cell lines and Sex and Gender in Research](#)

Cell line source(s)

Human multiple myeloma cell line MM.1S (ATCC #CRL-2974), chronic myelogenous leukemia cell line K562 (ATCC #CCL-243), Burkitt's lymphoma cell line Raji (ATCC #CCL-86), acute lymphoblastic leukemia cell line NALM-6 (ATCC #CRL-3273), melanoma cell line A375 (ATCC #CRL-1619), glioblastoma cell line T98G (ATCC #CRL-1690) and U87MG (ATCC #HTB-14),

	hepatocellular carcinoma cell line HEP3B (ATCC #HB-8064), and HEK 293T (ATCC #CRL-3216) were purchased from the American Type Culture Collection (ATCC). All cell lines used in this study were tested mycoplasma free.
Authentication	MM.1S, K562, Raji, NAML-6, A375, T98G, U87MG, HEP3B, HEK 293T cells were authenticated by ATCC. All the FG engineered cells and BCMA knockout cells were not authenticated.
Mycoplasma contamination	All cell lines were tested negative for mycoplasma contamination.
Commonly misidentified lines (See ICLAC register)	No misidentified lines were used.

Animals and other research organisms

Policy information about [studies involving animals](#); [ARRIVE guidelines](#) recommended for reporting animal research, and [Sex and Gender in Research](#)

Laboratory animals	NOD.Cg-PrkdcSCIDII2rgtm1Wjl/SzJ (NOD/SCID/IL-2R γ ^{-/-} , NSG) mice were maintained in the animal facilities of the University of California, Los Angeles (UCLA) under the following housing conditions: temperature ranging from 68°F to 79°F, humidity maintained at 30% to 70%, a light cycle of On at 6:00 am and Off at 6:00 pm, and room pressure set to negative. 6-10 weeks old female mice were used for all experiments unless otherwise indicated.
Wild animals	No wild animals were used.
Reporting on sex	6-10 weeks old female mice were used for all experiments unless otherwise indicated.
Field-collected samples	This study did not involve field-collected samples.
Ethics oversight	All animals were maintained at the UCLA animal facilities and all animal experiments were approved by the Institutional Animal Care and Use Committee of UCLA.

Note that full information on the approval of the study protocol must also be provided in the manuscript.

Flow Cytometry

Plots

Confirm that:

- The axis labels state the marker and fluorochrome used (e.g. CD4-FITC).
- The axis scales are clearly visible. Include numbers along axes only for bottom left plot of group (a 'group' is an analysis of identical markers).
- All plots are contour plots with outliers or pseudocolor plots.
- A numerical value for number of cells or percentage (with statistics) is provided.

Methodology

Sample preparation	Purified cord blood derived human CD34+ cells were purchased from the HemaCare. Healthy donor PBMCs were obtained from the UCLA/CFAR Virology Core Laboratory without identification information under federal and state regulations. Primary MM patient bone marrow samples were collected at the Ronald Reagan UCLA Medical Center from consented patients through an IRB-approved protocol (IRB#21-001444) and processed. AlloNKT and Allo/15CAR-NKT cells were generated by differentiating gene-engineered cord blood CD34+ HSCs in a 5-stage Ex Vivo HSC-Derived NKT Cell Culture. Healthy donor PBMCs were used to generate the PBMC-derived conventional $\alpha\beta$ T, NKT, and NK cells (denoted as PBMC-Tc, PBMC-NKT, and PBMC-NK cells, respectively). In mouse study, cells were collected from mouse organs (e.g., bone marrow, spleen, liver, lung, blood, and peritoneal cavity).
Instrument	MACSQuant Analyzer 10 flow cytometer (Miltenyi Biotech)
Software	FlowJo software version 9 (BD Biosciences)
Cell population abundance	For cell sorting, highest purity method was used in FACS Aria II Flow Cytometer and 98% post-sorting purity was verified through flow cytometry analysis.
Gating strategy	SSC-A/FSC-A, SSC-W/SSC-H and FSC-W/FSC-H were used to identify singlets. Fixable Viability Dye eFluor 506 (Biolegend) was used to exclude dead cells in flow analysis. DAPI (Thermo Fisher Scientific) was included to exclude dead cells in FACS sorting. Gating strategy for each cell subtype in each experiment was described in the article and provided as panels in supplementary figures.

- Tick this box to confirm that a figure exemplifying the gating strategy is provided in the Supplementary Information.